


### A & D TM-2430 24 Hour Ambulatory Blood Pressure Monitor



Actual Size of Monitor



- Free 'Doctor Pro Software' with analysis capabilities including Trend, Correlation, Histograms, Circadian Rhythm analysis and a Hyperbaric Index - HBI or BP load
- Compact lightweight 24-hour ambulatory blood pressure monitor
- 3 modes of measurement intervals: Estimated Preset Sleep Intervals, Actual Sleep Intervals and Variable Time Intervals
- Automatically adjusts: proper cuff pressure, exhaust & end of measurement
- Pulse correction algorithm for faster measurement & noise rejection
- Automatic or manual measurement
- Date and time display
- 350 measurement memory capacity
- Data can be sent to print or attached to printout files in EMIS, Torex etc.
- Dimensions: 10H x 7W x 2.7D cm
- Weight: 215g
- Supplied complete with adult cuff, large cuff, 2 cuff covers, 10 record cards, 6 x NiMH rechargeable batteries & charger and carry case



4532	TM-2430 24-Hour Ambulatory B.P. Monitor	£1250.00
4536	Replacement Large Cuff 28-36cm for TM-2430	£110.00
4537	Replacement Adult Cuff 20-31cm for TM-2430	£110.00

Welch Allyn ABPM 6100S 24 Hour  
Ambulatory Blood Pressure Monitor

WelchAllyn®



systemone  
ONEPATIENTONERECORD



Demo Available



Free Welch Allyn DS-55 Durashock on Demo request until 31st March 2011

Free Welch Allyn Diagnostic set with every purchase until 31st March 2011



Paperless ABPM data can attach to EMIS, Vision, Torex etc.

- Improved technology for enhanced 24 hour blood pressure monitoring
- Lightweight and quiet in use to improve patient compliance
- Easy to set up and operate
- Dynamic cuff inflation automatically limits cuff pressure for patient comfort
- 100% latex-free cuff to help reduce risk of allergic reaction
- Program through your PC to standard or customised measurement protocols
- Supplied with standard machine washable adult cuff
- Optional Cuff Set, code 8974, comprises of Small Adult, Adult and Adult Plus Cuffs



Latex Free

2 Year Warranty

BHS A/A Validation

Applies to Monitor only

6159	ABPM 6100S Monitor and Software	£1694.00
6160	ABPM 6100S Small Adult Cuff 21.1-26.6cm	£48.00
6161	ABPM 6100S Adult Cuff 25.5-34.3cm	£48.00
6162	ABPM 6100S Adult Plus Cuff 27-42cm	£62.00
6163	ABPM 6100S Large Adult Cuff 34.3-48.2cm	£75.00
8974	ABPM 6100S Cuff Set	£136.00



# blood pressure monitors

## CareFusion Pulse Trace PCA 2




- Complete pulse contour analysis system
- Designed for use in clinical practice, the PCA 2 is an effective tool for the early detection, treatment and ongoing management of hypertension and coronary arterial disease, targeting organ damage and endothelial dysfunction
- Uses accurate and noise-free digital volume pulse (DVP) waveform analyses using high fidelity photo-plethysmography with auto power adjustment
- The clear, touch-sensitive colour LCD display ensures that all parameters can be set easily and quickly
- The built-in thermal printer allows records of an examination to be produced immediately and held with a patient's records
- The unit can be connected to an external printer if a non-thermal paper record is required
- Convenient, high resolution touch screen
- 2000 Test memory
- Spot-check and automated three test protocol modes
- Battery or mains powered
- Supplied complete with all accessories in a sturdy case

9426 PulseTrace PCA 2

£2729.00

## BoSo Ankle Brachial Index Blood Pressure System ABI 100



Demo Available

Latex Free

Validated in Germany, quick and accurate measurement of differential blood pressure for easy diagnosis of impending stroke, heart attack, vascular abnormalities in diabetes, hypertension, smokers, obese and the elderly

- Allows PAOD to be identified with precision at an early stage
- Nurses and assistants can use the machine with minimal training
- Simplifies the procedure to such an extent that early diagnosis is possible even if symptoms are not yet experienced
- All results are stored in the software database for onward transmission or comparison over time
- Once measured the doctor is presented with other important figures such as individual blood pressure readings for all arms and legs as well as the left right pressure difference that can indicate subclavian artery stenosis. The device is also able to indicate the possible presence of Arrhythmia

7291	ABI System 100 with 4 cuffs and software	£2595.00
6871	Cliniline Trolley with 2 Shelves, Guard Rail and Drawer	£255.00

# blood pressure monitors

Omron 907 Automatic B.P. Monitor

**OMRON**



**Demo Available**



- Choice between oscillometric and auscultatory mode (mmHg)
- Pressure setting as optional function for individual use
- Sequential measure mode average readings display
- High speed inflation with optimum inflation level setting
- Adjustable deflation rate in auscultatory mode
- Check Mode for checking accuracy
- Latex free, hand washable cuff
- Durable cuff 60,000 measurements
- Silent device
- Rechargeable batteries up to 300 measurements between charges
- Supplied with mains adaptor, rechargeable battery, normal size adult cuff, 1m tube, instruction manual
- Weight: 910g

**Latex Free**

**3 Year Warranty**

**BHS A/A Validation**

4845	907 B.P. Monitor	£298.50
4850	907 Mobile Stand	£128.00
4846	907 Small Cuff 17-22cm	£30.00
4847	907 Medium Cuff 22-32cm	£39.00
4848	907 Large Cuff 32-42cm	£47.00
4853	907 Replacement 1m Tubing	£6.25
4854	907 Replacement Rechargeable Battery	£30.00



# blood pressure monitors

## Omron i-C10 B.P. Monitor



- Fully automatic upper arm blood pressure monitor
- Intellisense™ technology enables personalised inflation
- Pre-formed comfort cuff fits both medium and large arms ranging from 22 – 42 cm
- Memory capacity for 84 measurements with date and time for 2 users
- Shows if irregular pulses are detected and only displays reliable results
- Dual user switch
- Morning hypertension indication
- Large 3-fold display (Sys/Dia/Pulse)
- Averaging for morning, evening and auto measurements
- Clinically validated to international Protocol & BHS
- Supplied with comfort cuff in a built-in hard storage case

3 Year Warranty

Latex Free

7570	i-C10	£99.50
4832	Comfort Cuff 22-42cm	£20.00
4826	Mains Adaptor Positive	£19.00

## Omron MIT Elite Plus B.P. Monitor

**OMRON**



- Validated for use in pregnancy and pre-eclampsia
- Fully automatic upper arm blood pressure monitor
- Intellisense™ technology enables personalised inflation
- 90 memories & data/time indication, keep track of exact time of measurement
- Detects and displays irregular heartbeats as well as body movement
- Large 3-fold display (Sys/Dia/Pulse)
- PC connectivity, download facility
- Display backlight
- Clinically validated to international Protocol
- Supplied with medium cuff, instruction manual, USB cable, CD-ROM, batteries and soft case

3 Year Warranty

Latex Free

4843	MIT Elite Plus	£99.50
4830	MIT Elite Plus Large Cuff 32-42cm	£12.50
7530	Mains Adaptor Positive for MIT only	£19.00

## Omron 705CP II B.P. Monitor with Printer



- Fully automatic upper arm blood pressure monitor
- Intellisense™ technology enables personalised inflation
- Bio Information Link (BI) functionality
- Memory capacity for 28 readings
- Large 3-fold display (Sys/Dia/Pulse)
- Includes printer for hard copies
- Clinically validated to international Protocol & BHS
- Supplied with medium cuff, instruction manual, printer, cable, batteries and soft case

3 Year Warranty

Latex Free

4837	705CP II with Printer	£134.00
4825-SM	Small Cuff 17-22cm	£12.50
4825-ME	Medium Cuff 22-32cm	£12.50
4825-LG	Large Cuff 32-42cm	£12.50
4826	Mains Adaptor Positive	£19.00
4842	Printer Paper, pk/5 rolls	£8.50

## Omron R6 Wrist B.P. Monitor



- Fully automatic wrist blood pressure monitor
- Intellisense™ technology enables personalised inflation
- Cuff fits 13.5 – 21.5 cm wrists
- Memory capacity for 90 memories, applicable for 2 users
- Memory has date and time stamp
- Large display showing blood pressure and pulse
- Advanced Positioning Sensor ensures accurate positioning resulting in accurate measurements
- Clinically validated to International Protocol
- Supplied in storage case with batteries

3 Year Warranty

Latex Free

4844	R6 Wrist	£79.00
------	----------	--------

# blood pressure monitors

## Omron M6 B.P. Monitor

3 Year Warranty

Latex Free



- Fully automatic upper arm blood pressure monitor
- Instellisense™ technology enables personalised inflation
- 90 memory capacity with date and time stamp
- Irregular heartbeat indicator, also features a movement indicator
- Blood pressure level indicator, easily identifies possible hypertension
- Innovative dual check system for 100% accuracy, first sensor takes the measurement while the second sensor checks its accuracy
- Cuff wrapping guide for accurate wrapping, a symbol on the display will confirm if the cuff is wrapped correctly or incorrectly
- Large 3-fold display (Sys/Dia/Pulse)
- Clinically validated to international Protocol & BHS
- Supplied with medium cuff, batteries and soft case

4827	M6	£69.00
4825-SM	Small Cuff 17-22cm	£12.50
4825-ME	Medium Cuff 22-32cm	£12.50
4825-LG	Large Cuff 32-42cm	£12.50
4826	Mains Adaptor Positive	£19.00

## Omron M6 Comfort B.P. Monitor

OMRON

3 Year Warranty

Latex Free

NEW



- Fully automatic upper arm blood pressure monitor
- Instellisense™ technology enables personalised inflation
- 90 memory capacity with date and time stamp
- Irregular heartbeat indicator, also features a movement indicator
- Blood pressure level indicator, easily identifies possible hypertension
- Innovative dual check system for 100% accuracy, first sensor takes the measurement while the second sensor checks its accuracy
- Cuff wrapping guide for accurate wrapping, a symbol on the display will confirm if the cuff is wrapped correctly or incorrectly
- Large 3-fold display (Sys/Dia/Pulse)
- Clinically validated to international Protocol & BHS
- Supplied with comfort cuff, which fits medium and large arms

7571	M6 Comfort	£85.00
4832	Comfort Cuff 22-42cm	£20.00
4826	Mains Adaptor Positive	£19.00

## Omron M2 B.P. Monitor



3 Year Warranty

Latex Free

- Fully automatic upper arm blood pressure monitor
- Instellisense™ technology enables personalised inflation
- Memory stores the last 21 recorded readings
- Irregular heartbeat indicator
- Large 3-fold display (Sys/Dia/Pulse)
- Clinically validated to International Protocol and BHS
- Supplied with medium cuff, batteries and soft case

7569	M2	£44.50
4825-ME	Medium Cuff 22-32cm	£12.50
4825-LG	Large Cuff 32-42cm	£12.50
4826	Mains Adaptor Positive	£19.00

## Omron M7 B.P. Monitor



3 Year Warranty

Latex Free

- Fully automatic upper arm blood pressure monitor
- Instellisense™ technology enables personalised inflation
- Pre-formed comfort cuff fits both medium and large arms ranging from 22 to 42 cm
- Memory capacity for 90 memories
- Memory has date and time stamp
- Large 3-fold display (Sys/Dia/Pulse)
- Averaging mode
- Clinically validated (international Protocol & BHS)
- Supplied with comfort cuff, batteries and hard case

4831	M7	£89.00
4832	Comfort Cuff 22-42cm	£20.00
4826	Mains Adaptor Positive	£19.00



# blood pressure monitors

## A & D UA-621 Fully Automatic Budget B.P. Monitor



**NEW**

**IHB**

**Latex Free**

**5 Year Warranty**

- Automatic one-touch operation
- WHO Blood Pressure result classification indicator
- 30 Memories
- Average BP Result display
- Irregular Heartbeat feature (IHB)
- Slim-fit comfortable adult cuff with correct size indicator
- 2 Minute automatic switch off for battery saving
- Optional Mains Adaptor
- Supplied with Adult cuff, Soft Storage Case, 4 x AA batteries and full instructions

9429	UA-621	£34.00
4529	Adult Cuff 22-32cm	£12.50
4530	Large Cuff 32-45cm	£17.00
4531	Optional Mains Adaptor	£13.00

## A & D UA-853 & 851 B.P. Monitor



**ESH Validation**

**IHB**

**Latex Free**

**5 Year Warranty**

- New Timewise™ technology allowing easy reading of results for key morning period
- Large 90 reading memory with timed AM and PM dated and average readings (AM and PM is not available on UA-851)
- Clock display with 3 alarms (not available on UA-851)
- Smooth, fast measurement
- WHO Blood Pressure Classification indicator
- Irregular Heartbeat Indicator
- Optional Mains Adaptor
- Supplied with Adult Cuff, Cloth Cover, 4 x AA batteries and full instructions

1020	UA-853 with Clock	£68.00
4527	UA-851 without Clock	£59.50
4528	Small Cuff 18-22cm	£12.50
4529	Adult Cuff 22-32cm	£12.50
4530	Large Cuff 32-45cm	£17.00
4531	Optional Mains Adaptor	£13.00

## A & D UA-767 Plus B.P. Monitor



**BHS A/A Validation**

**IHB**

**Latex Free**

**5 Year Warranty**

- Automatic one-touch operation
- WHO Blood Pressure classification indicator
- Last reading memory
- Slim-fit comfortable cuff with correct size indicator
- Can not be used with small cuff
- Optional Mains Adaptor
- Supplied with Adult Cuff, Storage Case, 4 x AA batteries and full instructions

1021	UA-767 Plus	£59.50
4529	Adult Cuff 22-32cm	£12.50
4530	Large Cuff 32-45cm	£17.00
4531	Optional Mains Adaptor	£13.00

## A & D UA-787 Plus B.P. Monitor with Multi Cuff



**ESH Validation**

**IHB**

**Latex Free**

**5 Year Warranty**

- Multi fit cuff fitting both medium and large arms, range 22 to 43cm
- 60 memory – timed and dated - clock feature with alarm settings
- WHO blood pressure classification indicator
- Irregular Heartbeat feature
- Advanced fuzzy-logic for accuracy and optimum comfort – for faster and almost silent operation
- Optional Mains Adaptor
- Supplied with Multi fit Cuff Storage Case, 4 x AA batteries and full instructions

4526	UA-787 Plus	£68.00
4528	Small Cuff 18-22cm	£12.50
4531	Optional Mains Adaptor	£13.00

# blood pressure monitors

## A & D TM-2655P Fully Automatic Practice B.P. Monitor



Supplied  
with  
**FREE STAND**  
value £310

- Irregular Heart Beat IHB Indicator
- AAMI / AANSI SP10 Clinically Validated
- Ideal for QOF clinical points commitment
- Saves time and practice resources
- Waiting room location can reduce level of white coat hypertension
- Ergonomic design gives ideal position for accurate BP measurement
- Quick print out of results for patient files
- Optional AD-6121A Class III Medical Scale available

### AD-6121A Class III Medical Scale

- Connects to TM-2655P Waiting Room B.P. Monitor
- Prints Date, Time and Weight on TM-2655P
- Screen patient's blood pressure and weight at the same time
- Class III approved for professional medical use
- Weighs from 2 to 150kg
- Graduation 100g
- Dimensions: 31.4D x 31.5W x 5.7H cm
- Weight: 3kg

4534	TM-2655P Practice B.P. Monitor	£1995.00
4535	Printer Rolls, pk/5	£30.00
9435	AD-6121A Medical Scale	£395.00



## A & D UA-789XL Bariatric B.P. Monitor



- Extra Large Cuff, for 42 to 60cm arm circumference
- New high power pump for quick inflation
- Advanced ECEV for speed, comfort and for small cuffs
- 60 memory + average timed and dated memory
- Clock display and alarm feature
- WHO Guide to give users a guide to the result
- Illuminated start button
- Optional Mains Adaptor
- Supplied with Bariatric Cuff, Storage Case, AA batteries and full instructions

4525	UA-789XL Bariatric B.P. Monitor	£119.50
7395	Optional Mains Adapter for 789XL	£40.00

## A & D UM-101 Mercury-Free Sphygmomanometer



- Used in exactly the same as a conventional mercury sphygmomanometer
- Mercury-Free blood pressure monitor
- Auscultatory blood pressure measurement
- Digital readout for pulse rate
- Adjustable display column for easier analysis of results
- Battery powered using 2 x AA supplied batteries
- Supplied with Adult Cuff
- Weight 1kg

7525	UM-101 Desk Top	£110.00
7526	UM-101 Mobile	£250.00
7527-SM	Small Cuff 18-22cm	£35.00
7527-AD	Adult Cuff 22-32cm	£35.00
7527-LG	Large Cuff 32-45cm	£35.00



# blood pressure monitors

## Accoson Duplex Aneroid



Latex Free

- Made in UK and guaranteed accurate to BS EN 1060-1
- One-handed use, with thumb and fingertip control
- Lightweight diecast case
- Supplied in soft case with zip
- Range 0 to 300mm Hg
- Latex free inflation system
- Supplied with single tube adult size velcro cuff with coiled tubing
- Available with Blue or Black cuff - specify colour when ordering

3732 Duplex Aneroid - specify colour £42.70

## Accoson Limpet Combine Model



Latex Free

- Made in UK and guaranteed accurate to BS EN 1060-1
- Cuff with integral stethoscope for self use
- Range 0 to 300mm Hg
- Lightweight plastic moulded case
- Supplied in soft case with zip

4021 Limpet 1 Combination £45.00

## Accoson Limpet 1 Aneroid



Latex Free

- Made in UK and guaranteed accurate to BS EN 1060-1
- One-handed use suitable for both right and left hands
- Supplied in soft case with zip
- Range 0 to 300mm Hg
- Supplied with latex free adult velcro cuff with straight tubing and inflation bulb
- Limpet Cuff Set comprises Limpet 1 with 4 Velcro cuffs; 1 x Outsize Adult, 1 x Alternative Adult, 1 x Standard Adult, 1 x Child

3729 Limpet 1 Aneroid £30.75  
3731 Limpet 1 Aneroid with Cuff Set £67.95

## Accoson Limpet 2 Aneroid



Latex Free

- Made in UK and guaranteed accurate to BS EN 1060-1
- One-handed use for both right and left hands
- Supplied in soft case with zip
- Range 0 to 300mm Hg
- Supplied with latex free adult velcro cuff, coiled tubing and inflation bulb

3730 Limpet 2 Aneroid £32.25

# blood pressure monitors

## Accoson 6 inch Aneroid



- Made in UK and guaranteed accurate to BS EN 1060-1
- Large diameter white dial
- Wall model has swivel mounting bracket
- Desk model has convenient carrying handle
- Stand model with inclined gauge and adjustable height, with a stable base on five castors
- Cuff is stored behind the gauge
- Range 0 to 300mm Hg
- For additional cuffs use double tube cuffs on page 42
- Supplied with latex free inflation system with coiled tubing

3733	6 inch Desk	1	£47.70
3734	6 inch Wall	2	£44.00
3735	6 inch Mobile	3	£121.50

## Accoson Greenlight 300

Accuracy confidence self-calibration without mercury

*Recommended by the British Hypertension Society*

greenlight 300™

- Recommended manual device of the British Hypertension Society
- A manual device with self-calibration each time it is switched on giving accuracy better than +/-3mm Hg
- Moving light indicates cuff pressure with bright large scale display
- None of the fragility of an aneroid gauge
- Frequent calibration checks no longer required
- Cuff deflation rate displayed to aid good technique
- Range 0 to 300mm Hg
- Designed within the NHS to provide an accurate alternative to mercury
- Wall mounting and Stand kits available
- Latex free inflation system and coiled tubing
- Full details and video on [www.greenlight300.com](http://www.greenlight300.com)
- Supplied as desk model with adult velcro cuff
- For additional cuffs use double tube cuffs on page 42



3736	Desk Model with Adult Velcro Cuff	4	£146.25
3737	Wall Conversion Kit		£26.25
3738	Mobile Stand Conversion Kit	5	£78.00

Recommended by the British Hypertension Society  
 - visit [www.bhsoc.org](http://www.bhsoc.org) (Manual Device)  
 - visit [www.greenlight300.com](http://www.greenlight300.com) for additional information & video



# blood pressure monitors

## Accoson Dekamet Mercurial



- Made in UK and guaranteed accurate to BS EN 1060-1
- Range 0 to 300mm Hg
- Lightweight diecast case
- Shatterproof plastic manometer tube
- Mercury reservoir with isolation valve
- Steel hinge pivots
- Supplied with latex free adult velcro cuff with double tubing and inflation bulb



## Accoson Freestyle Mercurial



- Made in UK and guaranteed accurate to BS EN 1060-1
- Range 0 to 300mm Hg
- Shatterproof plastic manometer tube
- Mercury isolation valve fitted to reservoir
- Wall model has swivel mounting bracket
- Cuff is stored behind manometer
- Supplied with latex free adult velcro cuff coiled tubing and inflation bulb



3726 Dekamet Mercurial £55.00

3727 Freestyle Desk 1 £55.00  
3728 Freestyle Wall 2 £55.00

## Mercury Spills Kit



- Kit enables staff to deal safely and easily with mercury spills
- In addition re-usable Mercury Spills Warning Card, dampening spray, warning tapes and instructions are include
- Comprising of two Disposal Units each containing: Mercury aspirator (syringe), Scrapers, Amalgamating powder, Collection sponge, Protective apron, Pair of Nitrile Gloves

7295 Mercury Spills Kit £46.95  
7296 Mercury Spill Kit Replacement Disposal Unit, pk/4 £61.50

## Accoson Velcro Cuffs & Accessories



	Single Tube		Double Tube		Double Tube with bulb and valve	
	Code	Price	Code	Price	Code	Price
Standard Adult 9" x 5", 22.8 x 12.7cm	3753	£11.60	3746	£11.60	3739	£15.20
Alternative Adult 14" x 5", 35.4 x 12.7cm	3754	£16.20	3747	£16.20	3740	£19.80
Outsize Adult 15" x 7", 38.1 x 17.8cm	3755	£25.20	3748	£25.20	3741	£28.80
Child 7.5" x 4", 19.1 x 10.2cm	3756	£9.10	3749	£9.10	3742	£12.70
Young Child 6" x 3", 15.2 x 7.6cm	3757	£9.10	3750	£9.10	3743	£12.70
Infant 4" x 2", 10.2 x 5.1cm	3758	£9.10	3751	£9.10	3744	£12.70
Thigh 20" x 7", 51 x 17.8cm	3759	£28.40	3752	£28.40	3745	£32.00



### Other Accessories

3760 Inflation Bulb & Valve £4.15  
3761 Tube Connectors, pair £1.75  
3762 Long Coiled Tube for wall and stand models £7.80  
3763 Short Coiled Tube for desk models £4.70

# blood pressure monitors

## Welch Allyn 767 Aneroid Sphygmomanometer

WelchAllyn®

Latex Free

NEW

Lifetime Calibration Warranty



See page 45 for Cuff Range

- Combines the traditional benefits with improved aesthetics and a new contemporary design
- Laser engraved dial face provides unmatched accuracy
- Certified accurate to +/- 3 mm Hg
- Recessed dial face and high contrast pointer increase visibility and reduce parallax errors
- Jeweled movement contributes to a longer life
- Mercury free
- Built in storage basket on wall and desk models
- New 5 leg mobile stand provides stability and maneuverability on mobile model
- Supplied with adult FlexiPort® BP cuff that is antimicrobial treated to help prevent mildew and bacteria growth

9982	767 Desk Aneroid	£93.00
9980	767 Wall Aneroid	£133.00
9981	767 Mobile Aneroid	£252.00

## Welch Allyn Connex® ProBP™ 3400 Digital Blood Pressure Monitor

NEW

- A reliable alternative to manual sphygmomanometers, mercury columns, and home-grade digital blood pressure devices
- Small, handheld design with various mountable options makes this device perfect for your practice
- Compatible with blood pressure cuff sizes ranging from small child to thigh
- Lithium ion battery technology delivers 100 readings on a single charge with rapid recharge time
- Obtains accurate readings and minimises error codes, even in the presence of many cases of motion or weak pulse
- Mercury and latex free to ensure clinician and patient safety
- 50 reading internal memory
- SureBP models take the blood pressure reading on inflation and capture the reading in just 15 seconds, improving patient comfort



See page 45 for Cuff Range

9877	ProBP 3400 Handheld	£420.00
9878	ProBP 3400 Wall	£468.00
9879	ProBP 3400, SureBP Handheld	£517.00
9880	ProBP 3400, SureBP Wall	£565.00

Latex Free

BHS A/A Validation



# blood pressure monitors


## Welch Allyn DuraShock™ Blood Pressure Range

The complete range of Welch Allyn hand held aneroid sphygmomanometers now incorporate shock-resistant DuraShock Technology and come complete with Welch Allyn Flexiport cuffs.

With this technology you can count on Welch Allyn gauges to:

- Remain in calibration longer than traditional blood pressure gauges, providing you with more accurate readings
- Not need to be replaced as often as traditional aneroids, giving you a lower lifetime cost
- Capable of passing AAMI's shock resistant specification of surviving a 76cm drop onto a hard surface and remain in calibration

# WelchAllyn®



### Welch Allyn Durashock™ DS58 Platinum



**NEW**

**Lifetime Calibration Warranty**

Applies to Gauge only

**Latex Free**

- Super shock-resistant - can withstand up to 152cm drop and remain in calibration
- Jewel movement for long life
- Laser-engraved dial face for accuracy
- Premium inflation system for comfort
- Ergonomic inflation bulb enhances user comfort
- Supplied with adult FlexiPort™ BP Cuff and zipper case

9983 DS58 Platinum £119.00

### Welch Allyn DuraShock™ DS66 Gold



**15 year Warranty**

Applies to Gauge only

**Latex Free**

- Trigger style air release valve
- Laser engraved dial for better accuracy
- Jewel movement for long life
- Shock resistant movement
- Clinically validated to AAMI standards
- Ergonomic design for added comfort and easy operation
- Supplied with a standard adult Flexiport™ Cuff and soft case

6191 DS66 Gold £80.00

### Welch Allyn DuraShock™ DS56



**15 year Warranty**

Applies to Gauge only

**Latex Free**

- Thumbscrew air release valve
- Customisable name label for easy staff identification
- Laser engraved dial for better accuracy
- Shock resistant movement
- Clinically validated to AAMI standards
- Supplied with a standard adult Flexiport™ Cuff and soft case

6190 DuraShock™ DS56, Black £69.00

### Welch Allyn DuraShock™ DS65 Silver



**10 Year Warranty**

Applies to Gauge only

**Latex Free**


- Trigger style air release valve
- Laser engraved dial for better accuracy
- Shock resistant movement
- Clinically validated to AAMI standards
- Ergonomic design for added comfort and easy operation
- Supplied with a standard adult Flexiport™ Cuff and soft case

6192 DS65 Silver £67.00

# blood pressure monitors

## Welch Allyn DuraShock™ DS55

- Thumbscrew air release valve
- Customisable name label for easy staff identification
- Shock resistant movement
- Clinically validated to AAMI standards
- Supplied with a standard adult Flexiport™ Cuff and soft case



See below for Cuff Range

**Latex Free**

**10 Year Warranty**

Applies to Gauge only

6185-BK	DuraShock™ DS55, Black	£58.00
6185-BL	DuraShock™ DS55, Blue	£58.00
6185-RD	DuraShock™ DS55, Red	£58.00
6185-YE	DuraShock™ DS55, Yellow	£58.00

## Welch Allyn DuraShock™ DS54

**Latex Free**

- Thumbscrew air release valve
- Shock resistant movement
- One tube friction fit
- Clinically validated to AAMI standards
- Supplied with a standard adult Flexiport™ Cuff and soft case



See below for Cuff Range

**5 Year Warranty**

Applies to Gauge only

6184	DuraShock™ DS54, Black	£44.00
------	------------------------	--------

## Welch Allyn FlexiPort™ Tubeless BP Cuffs

- Easy one-hand attachment or detachment of tubing from FlexiPort™ makes changing cuffs much faster and easier
- Meets all the latest clinical guidelines for proper fit from AAMI and AHA
- Tested to withstand demanding clinical environments
- Colour-coded for easier visual size identification of cuffs

**Latex Free**

**3 Year Warranty**



6659	Child Cuff Green, 15–21 cm	£15.90
6660-SM	Small Adult Cuff Blue, 20–26cm	£19.50
6660-AD	Adult Cuff Black, 25–34cm	£19.50
6660-LG	Large Adult Cuff Burgundy, 32–43cm	£19.50
6665	Thigh Adult Cuff Black, 40–55cm	£23.00



# blood pressure monitors

## Guardian Pro Berlin



5 Year Warranty

Latex Free

Applies to Gauge only

- High quality aneroid with pressure safe mechanism
- Air release valve can be positioned to suit left or right handers
- Luminous dial for low light situations
- Impact resistant housing
- Chromium plated body
- Supplied with storage pouch and adult cuff

6946 Pro Berlin £38.00

## Guardian Pro Desk/Wall



5 Year Warranty

Latex Free

Applies to Gauge only

- Measurement range 0 to 30cm Hg
- Spiral Tube
- Accurate and reliable
- Easy to read large dial
- Supplied with adult cuff

6938 Desk Model 1 £64.50  
6939 Wall Model 2 £59.50

## Guardian Pro Mobile

GUARDIAN<sup>®</sup>  
EMERGENCY & RESUSCITATION



5 Year Warranty

Applies to Gauge only

Latex Free

- Measurement range 0 to 30cm Hg
- Spiral Tube
- Accurate and reliable
- Easy to read large dial
- Height adjustable non tip stand
- Supplied with storage basket and adult cuff

6940 Mobile Model £118.00

## Guardian Pro Europa



5 Year Warranty

Latex Free

Applies to Gauge only

- Swivel body, suitable for left and right handed operation
- High quality aneroid with pressure safe mechanism
- Air release valve can be positioned to suit left or right handers
- Easy read dial
- Impact resistant housing
- Chromium plated body
- Slip resistant inflation bulb
- Supplied with storage pouch and adult cuff

6947 Pro Europa £39.00

## Clip-on



5 Year Warranty

Latex Free

Applies to Gauge only

- Reliable sphygmo, ideal for nursing
- Separate dial and inflation bulb
- Dial clips to cuff for easy reading
- Screw valve air release
- Supplied with storage pouch and adult cuff

6948 Clip-on £22.00

## Clip-on Combination



5 Year Warranty

Latex Free

Applies to Gauge only

- Simple, economically priced self test BP device
- Individual dial and inflation systems
- Dial clips on to cuff for easy reading
- Stethoscope attached to cuff for instant arterial placement
- Supplied with storage pouch and adult cuff

6949 Clip-on Combination £25.50

# blood pressure monitors

## BoSo Classico



15 Year Warranty

Latex Free

Applies to Gauge only

- Metal construction in matt black finish
- Smooth air release system
- Twin tube inflation system, stops 'jumping' of needle
- Accurate and reliable
- Shock proof manometer
- Adult cuff
- Supplied in carry pouch

6941 Classico £51.95

## BoSo Varius



10 Year Warranty

Latex Free

Applies to Gauge only

- Economically Priced - With no compromise on quality
- Shock and Corrosion Protected Manometer
- Safe Against Overpressure - A common cause of manometer failure from over-inflation
- 6cm Dial - Clear and easy to read
- Unbreakable Plastic Coating
- Supplied in carry pouch

6942 Varius £33.00

## BoSo Profitest



10 Year Warranty

Latex Free

Applies to Gauge only

- Push button valve, instant air release
- Two-tone design
- Shockproof manometer
- Strong and light ABS construction
- Clear and easy to read 6cm dial
- Quick and easy to use
- Adult cuff
- Supplied in carry pouch

6943 Profitest £38.50

## BoSo Clinicus



10 Year Warranty

Latex Free

Applies to Gauge only

- Lightweight ABS construction
- Shockproof manometer
- Strong and durable
- Traditional air release system
- Large easy to read 6cm dial
- Supplied with carry pouch

0813 Clinicus £41.50

## Boso Roid



15 Year Warranty

Latex Free

Applies to Gauge only

- Precise easy read 6cm dial
- Thread lock cuff connector system
- Shockproof manometer
- Protective metal rim
- Traditional air release system
- Adult cuff
- Supplied in carry pouch

6945 Roid £49.70

## Boso Egotest



10 Year Warranty

Latex Free

Applies to Gauge only

- Combination aneroid sphygmomanometer/ stethoscope
- Automatic push-button air release
- Tough ABS construction
- Stethoscope with cuff fastening for instant arterial location
- Large 6cm dial
- Supplied with full user instructions and storage case

0815 Egotest £44.50



# blood pressure monitors

## Riester Big Ben Mobiles

**Riester**  
The familiar way



Latex  
Free



- Large easy-read dial
- Economically priced
- Specially tempered copper beryllium diaphragm
- Integral cuff storage basket
- Range 0 to 30cm Hg
- Accurate +/- 3mm Hg
- Supplied with adult cuff

5827	Big Ben Round Mobile	<b>1</b>	£137.00
5828	Big Ben Square Mobile	<b>2</b>	£137.00
5829	Big Ben Adult Cuff		£22.80
5830	Big Ben Child Cuff		£22.80
5831	Big Ben Obese Cuff		£39.90

## Riester Big Ben Desk/Wall Models

Latex  
Free



- Large easy-read dial
- Economically priced
- Accurate +/- 3mm Hg
- Specially tempered copper beryllium diaphragm
- Integral cuff storage basket
- Swivel facility for wall mounted version
- Range 0 to 30cm Hg
- Supplied with adult cuff

5823	Ben Round Desk Model	<b>3</b>	£74.00
5824	Big Ben Round Wall Model		£68.00
5825	Big Ben Square Desk Model		£74.00
5826	Big Ben Square Wall Model	<b>4</b>	£68.00
5829	Big Ben Adult Cuff		£22.80
5830	Big Ben Child Cuff		£22.80
5831	Big Ben Obese Cuff		£39.90

# blood pressure monitors

## Riester Shock-proof



- Shock proof from a height of 120cm
- Patent pending metal air channel technology for stable needle deflection
- Fast and ergonomic inflation using a unique design with integrated spoon
- Scale diameter 5.2cm White, with Blue needle
- Precise and wear-free air release valve with fine adjustment mechanism
- Metal screw connection for simple and fast change of cuff
- Maximum error tolerance +/- 3mm Hg
- Supplied with vinyl zipper case with adult cuff

5832 Shock-proof £41.00

## Riester R1 Shock-proof Set

**Riester**  
The familiar way



- Shock-proof up to a falling height of 120cm
- Scale: Ø 52mm white with blue pointer
- Valve base and air duct with metal precision valve
- Two-component ABS casing
- Metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomic operation
- Maximum error tolerance of +/- 3mm Hg
- Weight of R1 shock-proof: 173g
- Set includes storage box and 3 cuffs -large adults, adults, children and R1 shock-proof with white scale



9431 R1 Shock-proof set with 3 cuffs £160.00

## Riester Ri-San



- High quality palm-style
- Impact resistant plastic housing
- Large, easy to read aluminium scale of 65mm Dia and up to 300mm Hg
- Fingertip control of finely tuned air release valve
- Calibrated adult velcro cuff
- Maximum error tolerance of unit +/- 3mm Hg
- Supplied in zipper bag with adult cuff
- Available in Slate Grey, Saffron Yellow and Green, specify colour when ordering



5833 Ri-San - specify colour £35.00  
 5837 Ri-San Adult Cuff £22.80  
 5838 Ri-San Child Cuff £22.80  
 5839 Ri-San Obese Cuff £39.90

## Riester Babyphone®



- Bear imprint on the cuff
- Manometer with chromium-plated metal casing
- Stainless steel ergonomic spoon-shaped handle
- Metal manometer with aluminium scale 4.9cm diameter
- Latex free inflation bulb with chromium-plated precision air-release valve
- Supplied in padded nylon zipper case with 1 tube, 3 velcro cuffs; children, infants and new-born

5836 Babyphone® £79.50

