

Pelican PELIspec



**BUY 4+
packs at
£24.75
each**

- Individually packed
- Available in Virgin, Small, Medium, Medium Long and Broad, please specify size when ordering



5728 PELIspec, pk/25

£29.95

ComfiSpec® Disposable Vaginal Speculum



**BUY 4+
packs at
£19.50
each**

- Optimum 38 degree opening angle provides best access gynaecology procedures
- Textured and frosted handle provides superior grip and control
- Predictable fail safe feature is outside the patient
- Reduced danger of skin trapping
- Fast acting locking mechanism for reduced procedure durations and less patient trauma
- Professional medical grade plastic
- Available in Small, Medium, Medium-Long and Broad type, please specify size when ordering



9640 ComfiSpec® Vaginal Speculum, pk/25

£22.95

ProSpec Sims Disposable Plastic Vaginal Speculum



- Sterile individually packed single use Sims vaginal speculum
- Available in Small, Medium and Large, please specify size when ordering

Small	55 x 27/65 x 30 mm
Medium	65 x 30/70 x 32 mm
Large	75 x 32/80 x 36 mm



6878 Sims Plastic Vaginal Speculum, pk/20

£19.00

Sims Disposable Stainless Steel Vaginal Speculum



- Sterile individually packed single use Stainless Steel Sims vaginal speculum
- Available in Small, Medium and Large, please specify size when ordering



9557 Sims Stainless Steel Vaginal Speculum

£3.00

ProSpec Disposable Vaginal Speculum

ProSpec

- High quality single use cusco vaginal speculum
- Made from non-toxic plastic and with non-heavy metals avoiding brittleness
- Strong and simple to use
- Screw thread ratchet system, patient friendly
- Supplied sterile individually packed in boxes of 20
- Available in bulk packs for extra savings
- Available in Mini, Virgin, Small, Medium, Large and Medium Extra Long, please specify size when ordering



6754	Vaginal Speculum, pk/20		£19.00
SAVE MORE			
6877	Mini	pk/270	£240.00
6872	Virgin	pk/230	£204.00
6873	Small	pk/240	£213.00
6874	Medium	pk/180	£160.00
6875	Large	pk/160	£142.00
6876	Medium Extra Long	pk/170	£151.00

Cusco Disposable Vaginal Speculum

- Featuring longer beak
- Clearer plastic improving visibility
- Threaded nut with safety stop feature
- Thumb ridges on handles for maneuverability
- Available in Small, Medium and Large, please specify size when ordering



7450	Cusco Vaginal Speculum, pk/40		£19.95
------	-------------------------------	--	--------

Medispec Disposable Vaginal Speculum

- Sterile and individually wrapped
- Rounded ends for patient comfort
- Latex free
- Lockable
- Colour coded
- Available in Small (Blue), Medium (Red) and Large (Green), please specify size when ordering



9953	Medispec Vaginal Speculum, pk/100		£54.00
------	-----------------------------------	--	--------

Cusco Disposable Stainless Steel Speculum


- Stainless Steel
- Available in Small, Medium and Large, please specify size when ordering



6541	Cusco Speculum		£1.52
------	----------------	--	-------



ProSpec Disposable Hegar Dilators




- Plastic
- Choice of 5 sizes: Size 1-2, size 3-4, size 5-6, size 7-8 & size 9-10, please specify size when ordering

6761 Hegar Dilators, pk/5

£5.50

IUCD Removal Hook



- Stainless Steel

6533 IUCD Removal Hook £1.60

Cervix Dilators Double Ended 22cm



- Stainless Steel
- Available in 1/2, 3/4, 5/6, please specify size when ordering

9589 Cervix Dilators £1.28

Sound and Dilator Double Ended



- Plastic

9567 Sound and Dilator £1.55

ProSpec Sims Uterine Sound



- Plastic
- Curved with bulbous tip and graduated
- Graduated up to 24cm, total length 32cm

6881 Uterine Sound, pk/5 £5.50

Sound 25cm



- Plastic

9568 Sound 25cm £0.71

Coil Remover



- Plastic
- Angled tip with 3:2 teeth combination
- Tip has 'V' design to aid smooth coil retrieval

5737 Coil Remover £0.70

ProBrush Cervical Sterile Sampling Brush



- Long flexible stem and 'arrow' shaped head allow for faster, trouble-free sampling

3845 Cervical Brush, pk/25 £6.00

Gynobrush Special



- Ideal for gynaecological sample with a special designed head that tapers to a fine tip
- Approximately 20cm long and supplied clean but not sterile
- Not to be used for endoetrial sampling

3846 Gynobrush, pk/100 £32.50

ProSpec Thread Retriever



- 12 Micro Teeth - hook the thread effectively
- Dimensions - overall length 29.8cm - handle 20.4cm to toothed section 9.4cm

3847 Thread Retriever, pk/5 £10.25

Teales Vulsellum Tissue Forceps 25cm



- Plastic

9565 Teales Vulsellum £0.88


Teales Vulsellum Forceps



- Stainless Steel

6534	Forceps Straight 23cm	£3.40
6532	Forceps Curved 18cm	£3.00
6535	Forceps Curved 23cm	£3.40

Allis Tissue Forceps



- Stainless Steel

9563	Allis Tissue Forceps 15cm	£2.38
9564	Allis Tissue Forceps 23cm	£3.08

Sponge Holder 24cm



- Plastic

9512 Sponge Holder 24cm £0.65

Ramplays Sponge Holder



- Stainless Steel

6536	Sponge Holder 18cm	£2.64
6537	Sponge Holder 24cm	£3.60

Uterine Sims Scissors Curved 20cm



- Stainless Steel

9576 Uterine Sims Curved £2.95

Tenaculum Forceps



- Plastic
- Single toothed grip

5740 Tenaculum £1.25

Braun Tenaculum Forceps 24.5cm



- Stainless Steel

9561 Tenaculum £3.05

Scissors Plastic Handle



- Long handle specifically designed to cut coil threads
- Choice of 2 sizes

5739	Scissors 11.4cm	£0.41
9575	Scissors 22cm	£0.94

Mayo Scissors



- Stainless Steel

6519	Scissors Straight 14cm	£1.35
6520	Scissors Straight 19cm	£2.00
6521	Scissors Curved 19cm	£2.25

Dressing Scissors



- Stainless Steel
- Straight

6526	Scissors Sharp/Sharp	£0.83
0883	Scissors Sharp/Blunt	£0.83

Ayres Spatula



- Type S
- Wooden
- Clinically clean

3848 Ayres Spatula £3.55

IUCD Instrument Pack Plastic



- Single use pack comprising of
- 1 x Sponge Forceps
- 1 x Uterine Scissors
- 1 x Tenaculum Forceps
- 1 x Uterine Sound



9534 IUCD Instrument Set Plastic

£3.95

IUCD Instrument Pack




- Single use pack comprising of
- 1 x Mayo Scissors Curved 19cm
- 1 x Rampleys Sponge Holder 24cm
- 1 x Teales Vulsellum Curved 23cm
- 1 x Procedure Sheet/Sterile Field 45 x 50 cm



7517 IUCD Instrument Pack

£7.50

IUCD Instrument Pack with Galabin Sound



- Single use pack comprising of
- 1 x Sims Uterine Scissors
- 1 x Rampley Sponge Forceps
- 1 x Teales Vulsellum Forceps Curved
- 1 x Uterine Galabin Sound



9535 IUCD Instrument Pack with Galabin Sound

£8.55

IUCD Comprehensive Packs



- Choice of packs; with or without Uterine Sound

- Single use pack comprising of
- 1 x Tray
- 1 x 120ml Gallipot
- 1 x 700ml Kidney Dish
- 1 x Mayo Scissor Curved 19cm
- 1 x Rampley Sponge Holder 24cm
- 5 x Large Cotton Wool Balls
- 1 x Teales Vulsellum Curved 23cm
- 1 x Cosco Medium Speculum
- 1 x Drape 75 x 95cm
- 1 x Dressing Towel 38 x 40 cm
- 1 x Sterile Field
- 1 x Galabin Uterine Sound
- available with code 6546 only



6546 IUCD Pack with Sound

£13.95

6547 IUCD Pack without Sound

£12.95

IUCD Fitting Pack



- Single use pack comprising of
- 1 x Vaginal Speculum
- 1 x Uterine Sound
- 1 x Rampley Sponge Forceps
- 1 x Uterine Scissor Curved
- 1 x Teales Vulsellum Forceps
- 1 x Lubricating Gel Sachet 5g



9536 IUCD Fitting Pack

£10.20

IUCD Premium Fitting Pack




- Single use pack comprising of
- 1 x Tray
- 1 x Crepe Paper
- 1 x Towel Dressing
- 1 x Paper Tissue
- 1 x 120ml Gallipot
- 1 x Kidney Tray
- 5 x Large Cotton Wool Balls
- 1 x Rampley Forceps
- 1 x Uterine Sound
- 1 x Teales Vulsellum Forceps
- 1 x Vaginal Speculum



9537 IUCD Premium Fitting Pack

£12.25

Basic Vaginal Examination Pack



- Single use pack comprising of
 - 1 x 60ml Gallipot
 - 1 x 500ml Bowl
 - 8 x Large Cotton Wool Balls
 - 1 x Dressing Towel 43 x 38cm
 - 1 x Sterile Field



6515 Basic Vaginal Exam Pack £0.65

Vaginal Examination Pack



- Single use pack comprising of
 - 1 x 60ml Gallipot
 - 1 x 120ml Gallipot
 - 5 x X-Ray Det Swab
 - 1 x Large Maternity Pad
 - 1 x Fenestrated Drape 60 x 60 cm
 - 1 x Dressing Towel 38 x 40 cm
 - 1 x Sterile Field
 - 1 x Yellow Disposal Bag



6517 Vaginal Exam Pack £1.25

Obs & Gynae Pack



- Single use pack comprising of
 - 1 x Tray
 - 2 x 120ml Gallipot
 - 2 x 10 x10cm Non Woven Swab 4 ply
 - 1 x Large Maternity Pad
 - 5 x Large Cotton Wool Balls
 - 2 x Dressing Towel 38 x 40 cm
 - 1 x Sterile Field



6516 Obs & Gynae Pack £0.95

Maternity Pack



- Provides all the essential items needed to handle an emergency birth

- Sterile pack comprising of
 - 2 x Large Maternity Pad
 - 1 x Square Aborbency
 - 1 x Maternity Bag
 - 1 x Non Powdered Gloves
 - 1 x Baby Towelling
 - 1 x Maternity Pack Label
 - 1 x Umbilical Cord Scissors
 - 1 x Umbilical Cord Clamp

7538 Maternity Pack £24.00

HRT Pack



- Single use pack comprising of
 - 1 x Tray
 - 1 x Plastic Blue
 - 1 x Plastic Blue Tweezers
 - 1 x Scalpel No 11
 - 2 x Vinyl Gloves, Non Sterile, Powder Free, Medium
 - 1 x Waterproof Drape 50 x 50 cm with 6 x 6 cm hole
 - 5 x Woven Swabs, 10 x 10 cm, 8 ply
 - 2 x Dressing Towel 43 x 38 cm
 - 1 x Trochar and Cannula
 - 1 x Plaster
 - 1 x Tisept
 - 1 x Injection Wipes
 - 1 x 5ml Syringe
 - 1 x 23G Hypodermic Needle
 - 1 x 21G Hypodermic Needle
 - 5 x Small Cotton Wool Balls

9538 HRT Pack £17.95

Implant Insertion Pack



- Single use pack comprising of
 - 1 x Blister With 2 Integral Gallipot
 - 1 x Fold-out Sterile Field 60 x 60 cm
 - 1 x Adhesive Fenestrated Sheet with Hole
 - 1 x Alcohol Pad 70
 - 1 x Syringe 2ml
 - 1 x Hypodermic Needle 16 x 0.5
 - 1 x Hypodermic Needle 25 X 0.8
 - 1 x Stretch Bandage
 - 2 x Non Woven Swabs 10 x 10 cm 4 ply
 - 1 x Plaster 5 x 7.2 cm
 - 1 x Yellow Bag & Tie



9539 Implant Insertion Pack £2.25

Implant Removal Pack



- Single use pack comprising of
 - 1 x Blister with 2 Integral Gallipot
 - 1 x Fold-out Sterile Field 60 x 60 cm
 - 1 x Adhesive Fenestrated Drape with Hole
 - 1 x Halstead Mosquito Forceps Curved
 - 1 x Scalpel No 15
 - 2 x Non Woven Swabs 10 x 10 cm 4 ply
 - 3 x Reinforced Skin Closures
 - 1 x Stretched Bandage
 - 1 x Syringe 2ml
 - 1 x Hypodermic Needle 16 x 0.5
 - 1 x Yellow Bag & Tie

9540 Implant Removal Pack £4.45

Cusco Vaginal Specula



- Available in Virgin, Small, Medium, Medium-Long, Large - please specify size when ordering

6492 Cusco Vaginal Specula £8.00

Winterton Vaginal Specula



6497 Small Long 13 x 2.6cm £24.50

Sims Vaginal Specula



- Double ended

6498 Medium 7 x 3.5cm £7.65

Hegar Dilators




- Individual sizes available
Size 1/2
3/4
5/6

6488 Set of 8, 1 x size - 3/4 to 17/18 £49.00

6487 Individual size - please specify £6.50

Braun Tenaculum Forceps



6489 Braun 25cm £11.65

Teale Vulsellum Forceps



6490 Teale 3 x 4 Teeth 23cm £8.25

IUCD Hook



6491 Removal Hook £4.10

IUCD Instrument Set



One each of:

- Sims Uterine Scissors COF 8"
- Spencer Wells Artery Forceps Str 8"
- Rampley Sponge Forceps 9 1/2"
- Luer Vulsellum Forceps 1 x 1
- Medium Cusco Vaginal Speculum
- IUCD Removal Hook
- Galabin Uterine Sound

6499 IUCD Set

£39.90

Sims Uterine Scissors



6451 Sims 20cm 8" £4.95

Check-Mate Pregnancy Test



- One-step immunochromatographic assay for the determination of hCG in urine
- Single use Test strips are sealed in
- Accuracy of more than 99% found testing random urine samples, compared to commercially available hCG IRMA (Serono MAIA) kit
- Built-in control is included in the test
- Dip-strips can be used to determine pregnancy as early as the first day of the missed period

7388 Pregnancy Test, pk/20 £24.00

Clearview Easy hCG



- Test for the detection of hCG in urine as an aid to the diagnosis of pregnancy
- Reliable and accurate results in 3 minutes
- Sensitive to 25mIU/ml
- Minimal training required
- Flexible as four sampling options possible
- Can detect pregnancy on day of missed period
- Internal control indicates if test has functioned correctly
- Room temperature storage so no need to refrigerate

9505 Clearview Easy hCG, pk/20 £27.95

K-Y Jelly



- Water soluble lubricating jelly

7543 K-Y Jelly 42g tube £2.30
7544 K-Y Jelly 82g tube £3.75

ComfiTest™ Pregnancy Test



- Early detection of pregnancy using either urine or serum
- Detects the presence of hCG at the sensitivity of 25mIU/mL
- Clinical studies show > 99% accuracy
- The test is capable of detecting pregnancy as early as 1 day after the first period missed
- Quick and easy to use, results can be read after only 3 minutes for urine and 5 minutes for serum
- Only 3 drops of urine or serum required to conduct the test

SAVE MORE	1	4+
9504 ComfiTest Pregnancy Test, pk/20	£24.95	£20.95

Instillagel® Anaesthetic Lubricant



Instillagel

- Sterile gel containing a local anaesthetic and antiseptic
- Supplied in syringe ready for use
- **Pharmacy item, can only be supplied to registered medical practitioner**

Instillaquill

- Sterile single use Instillaquill Applicator Extension Tubes facilitate the application of Instillagel for virtually any gynaecological procedure

7561 Instillagel® Anaesthetic Lubricant 6ml, pk/10	1	£16.50
8167 Instillagel® Anaesthetic Lubricant 11ml, pk/10		£18.75
7562 Instillaquill, pk/10	2	£32.00

ComfiGel®



- May be used where a sterile, water soluble, non-staining lubricating jelly is indicated
- Non-irritating to the skin, tissues and mucus membranes
- Does not affect surgical instruments, rubber or plastics
- Available in 42g tube, 82g tube and boxes of 100 x 5g sachets

SAVE MORE	1	4+
9512 ComfiGel®, 42g tube	£1.55	£1.35
9513 ComfiGel®, 82g tube	£1.95	£1.69
9514 ComfiGel®, 5g sachets, pk/100	£21.25	£19.50

