

# ophthalmoscopes & otoscopes

## Welch Allyn 2.5v PocketScope™ Set

**WelchAllyn®**

**10 Year Warranty**



### 2.5v Pocket Ophthalmoscope Head:

- 6 diagnostic beams: Small, Medium and Large Spot, Target, Slit, Cobalt Blue and Red-free Filter
- Dioptre range: +40 to -25
- Double sealed, dust-free optics
- Fibre optic otoscope with 2.5 x magnification

### 2.5v Pocket Otoscope Head:

- Fibre optic 2.5x magnification
- Insufflation port for pneumatic otoscopy
- Supplied with 1 x AA metal handle, pocket clip, 4 reusable specula and 20 disposable specula in a hard case - batteries not supplied

### Ophthalmoscope Lenses




|      |   |         |
|------|---|---------|
| 6138 | Pocketscope™ Set  | £213.00 |
| 6805 | Duracell AA Size Batteries, pk/4                                    | £2.90   |
| 6645 | Bulb 03900 for 2.5v Pocket Ophthalmoscope                           | £14.90  |
| 6646 | Bulb 03400 for 2.5v Pocket Otoscope                                 | £14.90  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 | £52.00  |

## Welch Allyn 3.5v Prestige Set



**1**

### Ophthalmoscope Lenses



Nearly full view of the Tympanic Membrane with Macroview

### 3.5v Elite Ophthalmoscope Head:

- Unique coaxial ophthalmoscope allows easy entry to the eye and shadow-free viewing of the retina
- Double sealed, dust-free optics
- Dioptre range: +40 to -25

### 3.5v Macroview Otoscope Head:

- MacroView otoscope has 4.2x magnification - 30% greater than a traditional otoscope
- Twice the field of view - see almost all the tympanic membrane without panning
- Adjustable focusing - suits different ear canal lengths and compensates for far sighted eyes
- Insufflation port for pneumatic otoscopy
- Tip grip for secure fastening and ejection of ear specula

- Supplied with C-Cell battery handle or Lithium Ion battery handle and 4 reusable specula in a hard case

**LED Bulb Option**

See Bulb Offer on page 192

**10 Year Warranty**

|      |   |                  |
|------|---|------------------|
| 6141 | Prestige Set with C-Cell handle                                     | <b>1</b> £272.00 |
| 8975 | Prestige Set with Lithium Ion handle                                | £553.00          |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope                            | £18.50           |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope                        | £72.00           |
| 6644 | Bulb 06500 for 3.5v MacroView Otoscope                              | £16.40           |
| 9695 | LED Bulb 06500 for 3.5v MacroView Otoscope                          | £64.00           |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 | £52.00           |

# ophthalmoscopes & otoscopes

## Welch Allyn 3.5v Professional Diagnostic Set

WelchAllyn®



### 3.5v Standard Ophthalmoscope Head:

- 6 diagnostic beams: Small, Medium and Large Spot, Target, Slit and Red-free Filter
- Dioptre range: +40 to -25

### 3.5v Standard Otoscope Head:

- Fibre optic with HPX bulb for truer tissue colour
  - 2.5x magnification
  - Insufflation port for pneumatic otoscopy
- Supplied with C-Cell battery handle, batteries and 4 reusable specula in a hard case



|      |   |         |
|------|---|---------|
| 6139 | Professional Set  | £187.00 |
| 6649 | Bulb 03000 for 3.5v Standard Ophthalmoscope                         | £14.90  |
| 6648 | Bulb 03100 for 3.5v Standard Otoscope                               | £15.90  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 | £52.00  |

Ophthalmoscope Lenses



## Welch Allyn 3.5v Elite Diagnostic Set



See Bulb Offer on page 192

### 3.5v Elite Ophthalmoscope Head:

- Unique coaxial ophthalmoscope allows easy entry to the eye and shadow-free viewing of the retina
- 6 diagnostic beams: small, medium & large spot, target, slit and cobalt blue
- Red-free, and polarising filters
- Dioptre range: +40 to -25
- Double sealed, dust-free optics

### 3.5v Elite Otoscope Head:

- Fibre optic otoscope with HPX bulb for truer tissue colour
  - 2.5x magnification and integral throat illuminator
  - Insufflation port for pneumatic otoscopy
- Supplied with C-Cell battery handle or Lithium Ion battery handle and 4 reusable specula in a hard case



|      |   |                  |
|------|---|------------------|
| 6140 | Elite Set with C-Cell handle  | <b>1</b> £225.00 |
| 8976 | Elite Set with Lithium Ion handle & transformer                     | £465.00          |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope                            | £18.50           |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope                        | £72.00           |
| 6648 | Bulb 03100 for 3.5v Elite Otoscope                                  | £15.90           |
| 9694 | LED Bulb 03100 for 3.5v Elite Otoscope                              | £62.00           |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 | £52.00           |

Ophthalmoscope Lenses



# ophthalmoscopes & otoscopes

## Welch Allyn 3.5v Rechargeable Elite Desk Set

**WelchAllyn®**



See Bulb Offer on page 192



1

Ophthalmoscope Lenses



### Desk Unit Charger:

- With 2 x NICAD rechargeable handles in a charger



2

### 3.5v Elite Ophthalmoscope Head:

- Unique coaxial ophthalmoscope allows easy entry to the eye and shadow-free viewing of the retina
- 6 diagnostic beams: Small, Medium & Large Spot, Target, Slit and Cobalt Blue
- Red-free and polarising filters
- Dioptre range: +40 to -25
- Double sealed, dust-free optics

### 3.5v Elite Otoscope Head:

- Fibre optic otoscope with HPX bulb for truer tissue colour
- 2.5x magnification and integral throat illuminator

### Elite Desk with Lithium Handles:

- Same ophthalmoscope and otoscope as above but with Lithium Ion handles
- Two Lithium handles in charger offer a lightweight, longer lasting 'on time' alternative to NICAD handles with no 'memory' effect

|      |   |   |         |
|------|---|---|---------|
| 6135 | Elite Desk Set  | 1 | £683.00 |
| 6632 | Elite Desk Set with Lithium Ion Handles                             | 2 | £952.00 |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope                            |   | £18.50  |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope                        |   | £72.00  |
| 6648 | Bulb 03100 for 3.5v Elite Otoscope                                  |   | £15.90  |
| 9694 | LED Bulb 03100 for 3.5v Elite Otoscope                              |   | £62.00  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 |   | £52.00  |

## Welch Allyn 3.5v Rechargeable Prestige Desk Set



See Bulb Offer on page 192



3

Ophthalmoscope Lenses



### Desk Unit Charger:

- With 2 x NICAD rechargeable handles in a charger



4

### 3.5v Elite Ophthalmoscope Head:

- Unique coaxial ophthalmoscope allows easy entry to the eye and shadow-free viewing of the retina
- 6 diagnostic beams: Small, Medium & Large Spot, Target, Slit and Cobalt Blue
- Red-free and polarising filters
- Dioptre range +40 to -25
- Double sealed, dust free optics

### 3.5v Macroview Otoscope Head:

- 4.2x magnification
- Twice the field of view - see almost all the tympanic membrane without panning
- Adjustable focusing
- Insufflation point for pneumatic otoscopy
- Tip grip firmly holds specula in position and ejection mechanism for hands-free disposal
- Throat illuminator

### Prestige Desk Set with Lithium Handles:

- Same ophthalmoscope and MacroView otoscope as above but with Lithium Ion handles
- Two Lithium handles in charger offer a lightweight, longer lasting 'on time' alternative to NICAD handles with no 'memory' effect

|      |   |   |          |
|------|---|---|----------|
| 6134 | Prestige Desk Set   | 3 | £750.00  |
| 6633 | Prestige Desk Set with Lithium Ion Handles                          | 4 | £1012.00 |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope                            |   | £18.50   |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope                        |   | £72.00   |
| 6644 | Bulb 06500 for 3.5v MacroView Otoscope                              |   | £16.40   |
| 9695 | LED Bulb 06500 for 3.5v MacroView Otoscope                          |   | £64.00   |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850 |   | £52.00   |

# ophthalmoscopes & otoscopes

## Welch Allyn 3.5v Elite Wall Unit

WelchAllyn®



See Bulb Offer on page 192

### Wall Unit:

- With 12 foot cables
- Handle based rheostats allow adjustment of light intensity while in use
- OptiSense™ automatic sensors turn handles on when removed and off when replaced

### 3.5v Elite Ophthalmoscope Head:

- Unique coaxial ophthalmoscope allows easy entry to the eye and shadow-free viewing of the retina
- Dioptre range: +40 to -25
- Double sealed, dust-free optics

### 3.5v Elite Oscope Head:

- Fibre optic otoscope with HPX bulb for truer tissue colour
- 2.5x magnification and integral throat illuminator



|      |   |         |
|------|---|---------|
| 6136 | Elite Wall Unit   | £612.00 |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope                          | £18.50  |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope                      | £72.00  |
| 6648 | Bulb 03100 for 3.5v Elite Oscope                                  | £15.90  |
| 9694 | LED Bulb 03100 for 3.5v Elite Oscope                              | £62.00  |
| 6651 | Disposable Oscope specula 2.75mm or 4.25mm - specify size, pk/850 | £52.00  |

Ophthalmoscope Lenses



## Welch Allyn 3.5v 767 Integrated Diagnostic Unit

A convenient, modular package that puts everything needed for the fast, accurate diagnosis at arm's length. By combining a wall aneroid, ophthalmoscope, otoscope, specula dispenser and thermometer (optional extra), Welch Allyn puts all the basic diagnostic tools within easy reach.

### Integrated 767 Diagnostic Wall Unit:

- Wall Transformer with automatic on/off optical sensors
- Optional locking collars are available to provide security for instrument heads
- Disposable specula dispenser with 1,000 specula
- TycoS Aneroid Sphygmomanometer is certified accurate to within 1%, with a lifetime calibration and manufacturing warranty on manometer
- Easy to install and comes on a sturdy predrilled mounting board

### 3.5v Elite Ophthalmoscope Head:

- 6 apertures (small, medium and large spot, target, slit and cobalt blue) Polarising and red-free filters plus 28 lenses give diagnostic flexibility

### 3.5v Elite Fibre Optic Oscope Head:

- Insufflation port for pneumatic otoscopy
- HPX bulb gives 30% more light than traditional halogen, showing truer tissue colour



See Bulb Offer on page 192

|      |  |   |         |
|------|--|---|---------|
| 6137 | 767 Integrated Diagnostic Unit               | 1 | £952.00 |
| 6179 | Suretemp Plus 692 Wall Mounted Thermometer   | 2 | £236.00 |
| 6182 | Probe Covers, pack of 1000                   |   | £35.00  |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope     |   | £18.50  |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope |   | £72.00  |
| 6648 | Bulb 03100 for 3.5v Elite Oscope             |   | £15.90  |
| 9694 | LED Bulb 03100 for 3.5v Elite Oscope         |   | £62.00  |

Ophthalmoscope Lenses



# ophthalmoscopes & otoscopes

Welch Allyn LED Bulbs for Ophthalmoscopes and Otoscopes

**WelchAllyn®**

**NEW**

**7 Year Warranty**

## LED Special Offer

*Purchase an LED bulb with your diagnostic set and we will refund you the cost of the halogen bulbs(s) that come with your purchase*

### 30X Longer Life

- Tested to 30X halogen lamps



### Energy Efficient

- 3X the on-time versus halogen with Welch Allyn 3.5V rechargeable handles
- Lower operating temperatures intended for greater patient/clinician comfort



- **Long life** backed by a 7-year warranty ensuring peace of mind and virtually eliminating inconvenient interruptions for lamp replacement
- **SureColor™ LED technology** enabling reduction of intensity to alleviate glare or enhance patient comfort while preserving the colour and clarity of the view
- **High energy efficiency** resulting in 3X the examination time with each charge of a Welch Allyn rechargeable handle



|      |  |          |        |
|------|--|----------|--------|
| 9695 | LED Bulb 06500 for MacroView Otoscope        | <b>1</b> | £64.00 |
| 9694 | LED Bulb 03100 for 3.5v Elite & Pro Otoscope | <b>2</b> | £62.00 |
| 9696 | LED Bulb 04900 for 3.5v Elite Ophthalmoscope | <b>3</b> | £72.00 |
| 9697 | LED Bulb 03800 for PanOptic Ophthalmoscope   | <b>4</b> | £80.00 |

## Welch Allyn Bulbs for Ophthalmoscopes and Otoscopes plus Consumables



### Replacement Bulbs

|      |   |          |        |
|------|---|----------|--------|
| 6644 | Bulb 06500 for MacroView Otoscope                   | <b>5</b> | £16.40 |
| 6645 | Bulb 03900 for Pocket Ophthalmoscope                |          | £14.90 |
| 6646 | Bulb 03400 for Klinik Plus & Pocket Otoscope        |          | £14.90 |
| 6647 | Bulb 04900 for 3.5v Elite Ophthalmoscope            | <b>6</b> | £18.50 |
| 6648 | Bulb 03100 for 3.5v Elite & Pro Otoscope            | <b>7</b> | £15.90 |
| 6649 | Bulb 03000 for 3.5v Halogen Standard Ophthalmoscope | <b>8</b> | £14.90 |
| 6650 | Bulb 03800 for PanOptic Ophthalmoscope              | <b>9</b> | £21.00 |

### Otoscope Speculum

|      |   |           |        |
|------|---|-----------|--------|
| 6640 | Set of 4 reusable specula to fit all WelchAllyn Otoscopes Sizes 2, 3, 4 & 5mm | <b>10</b> | £15.00 |
| 6635 | Individual Size Reusable specula 2, 3, 4, or 5mm - specify size               | <b>11</b> | £5.10  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm - specify size, pk/850           | <b>12</b> | £52.00 |

|      |   |           |         |
|------|---|-----------|---------|
| 6643 | 3.5v Nicad Rechargeable Battery Orange, 72300 | <b>13</b> | £64.00  |
| 6642 | 3.5v Nicad Rechargeable Battery Black, 72200  | <b>14</b> | £60.00  |
| 6641 | 3.5v Lithium Ion Rechargeable Battery, 71960  | <b>15</b> | £127.00 |

# ophthalmoscopes & otoscopes

## Welch Allyn 2.5v PocketScope™ Otoscope

WelchAllyn®

- Fibre optic otoscope
- 2.5x magnification
- Insufflation port for pneumatic otoscopy
- Supplied with AA metal handle, pocket clip, throat illuminator, 4 reusable specula and 20 disposable specula in a hard case - batteries not supplied

10 Year Warranty



|      |  |         |
|------|--|---------|
| 6631 | PocketScope Otoscope   | £122.00 |
| 6805 | Duracell AA Size Batteries, pack of 4                                  | £2.90   |
| 6646 | Bulb 03400 for 2.5v Pocket Otoscope                                    | £14.90  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm<br>- specify size, pk/850 | £52.00  |

## Welch Allyn 3.5v Professional Otoscope

LED Bulb Option

See Bulb Offer on page 192

- Fibre optic otoscope with HPX™ bulb for truer tissue colour
- 2.5x magnification
- Insufflation port for pneumatic otoscopy
- Supplied with C-Cell battery handle, batteries and 4 reusable specula in a hard case

10 Year Warranty



|      |  |         |
|------|--|---------|
| 6634 | Professional Otoscope  | £128.00 |
| 6648 | Bulb 03100 for 3.5v Professional Otoscope                              | £15.90  |
| 9694 | LED Bulb 03100LED for 3.5v Elite Otoscope                              | £62.00  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm<br>- specify size, pk/850 | £52.00  |

## Welch Allyn 3.5v MacroView™ Otoscope

Nearly full view of the Tympanic Membrane with Macroview



- 4.2x magnification – 30% greater than a traditional otoscope
- Twice the field of view – see almost all the tympanic membrane without panning
- Adjustable focusing – suits different ear canal lengths and compensates for far-sighted eyes
- Insufflation port for pneumatic otoscopy
- Tip grip for secure fastening and ejection of ear specula
- Supplied with C-Cell battery handle, batteries and 4 reusable specula in a hard case

LED Bulb Option

See Bulb Offer on page 192

10 Year Warranty



|      |  |         |
|------|--|---------|
| 6133 | MacroView Otoscope   | £190.00 |
| 6644 | Bulb 06500 for 3.5v MacroView Otoscope                                 | £16.40  |
| 9695 | LED Bulb 06500 for MacroView Otoscope                                  | £64.00  |
| 6651 | Disposable Otoscope specula 2.75mm or 4.25mm<br>- specify size, pk/850 | £52.00  |



# ophthalmoscopes & otoscopes

## Keeler 2.8v Jazz Ophthalmoscope



**Keeler**



### Ophthalmoscope Lenses



- Brilliant Xenon illumination
- +/- 20D range of powers
- Renowned high performance Keeler optics with aspherical lens
- Total of 6 diagnostic beams
- Supplied with 3 interchangeable colour rings in a hard case
- Requires 2 x AA batteries - not supplied

|      |                                    |         |
|------|------------------------------------|---------|
| 9956 | Jazz Ophthalmoscope 2.8v           | £120.00 |
| 6805 | Duracell AA Size Batteries, pk/4   | £2.90   |
| 5806 | Bulb for Jazz Ophthalmoscope, pk/2 | £50.00  |

## Keeler 2.8v/3.6v Standard Ophthalmoscope



### Ophthalmoscope Lenses



### Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

- The positive maywheel lens control allows you to click through lens powers from +40D to -25D allowing detailed increased lens range allows detailed examination of the ocular structures
- Two diagnostic beams - a wide angle beam for general examination and a macular beam for more detailed examination of the central retina
- Red free filter for enhanced vessel examination can be used with both diagnostic beams, particularly useful for diabetic assessment and monitoring
- The integral soft brow rest allows for easy positioning and protects your lenses if you are a spectacle wearer
- Supplied in hard case with spare bulb
- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied

|      |   |         |
|------|---|---------|
| 5766 | Standard Ophthalmoscope 2.8v              | £179.00 |
| 5767 | Standard Ophthalmoscope 3.6v Rechargeable | £264.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 5797 | Lithium Mini Charger                      | £97.00  |
| 5807 | Bulb for 2.8v Ophthalmoscope, pk/2        | £50.00  |
| 5808 | Bulb for 3.6v Ophthalmoscope, pk/2        | £50.00  |

## Keeler 2.8v/3.6v Practitioner Ophthalmoscope



### Ophthalmoscope Lenses



### Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

- The positive maywheel lens control allows you to select from a range of lens powers from +40D to -25D. This allows increased lens range allows detailed examination of the ocular structures and guarantees ease of use for spectacle wearers
- The Practitioner includes 6 beams for a full diagnostic capability and a blue filter, particularly useful for assessing corneal scratches in children
- Red free filter for enhanced vessel examination can be used with both diagnostic beams - particularly useful for diabetic assessment and monitoring
- Supplied in hard case with spare bulb
- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied

|      |   |         |
|------|---|---------|
| 5768 | Practitioner Ophthalmoscope 2.8v              | £197.00 |
| 5769 | Practitioner Ophthalmoscope 3.6v Rechargeable | £281.50 |
| 6805 | Duracell AA Size Batteries, pk/4              | £2.90   |
| 5797 | Lithium Mini Charger                          | £97.00  |
| 5807 | Bulb for 2.8v Ophthalmoscope, pk/2            | £50.00  |
| 5808 | Bulb for 3.6v Ophthalmoscope, pk/2            | £50.00  |

# ophthalmoscopes & otoscopes

## Keeler 2.8v Jazz Fibre Optic Otoscope



**Keeler**



**NEW**

- Brilliant Xenon illumination
- 3x magnifying lens
- Supplied with 3 interchangeable colour rings in a hard case
- Requires 2 x AA batteries – not supplied
- Disposable specula available in 2.5, 3.5, 4.5 and 5.5mm

|      |   |        |
|------|---|--------|
| 9957 | Jazz F.O. Otoscope 2.8v                   | £85.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90  |
| 5811 | Bulb for Jazz Otoscope, pk/2              | £50.00 |
| 5791 | Disposable specula, pk/500 – specify size | £49.00 |

## Keeler 2.8v/3.6v Standard Otoscope



Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

- Closed head design particularly suited to pneumatic testing
- 2.5x magnifying lens removable to allow the introduction of small instruments
- Supplied in a hard case with 4 x reusable specula spare bulb, pneumatic testing adaptor and bulb extractor
- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

|      |  |         |
|------|--|---------|
| 5771 | Standard Otoscope 2.8v                       | £125.50 |
| 5772 | Standard Otoscope 3.6v Rechargeable          | £210.50 |
| 6805 | Duracell AA Size Batteries, pk/4             | £2.90   |
| 5797 | Lithium Mini Charger                         | £97.00  |
| 5811 | Bulb for 2.8v Pocket Standard Otoscope, pk/2 | £50.00  |
| 5812 | Bulb for 3.6v Pocket Standard Otoscope, pk/2 | £50.00  |
| 5791 | Disposable specula, pk/500 – specify size    | £49.00  |

## Keeler 2.8v/3.6v Practitioner Otoscope



Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

- Unobstructed open style head for easy introduction of small instruments for wax or foreign body removal
- Large slide aside 2.5x lens for clear view of aural structures and procedures under magnification
- Facility to undertake pneumatic testing for tympanic mobility
- Supplied in a hard case with 4 x reusable specula spare bulb and airtight chamber/pneumatic testing adaptor
- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

|      |  |         |
|------|--|---------|
| 5773 | Practitioner Otoscope 2.8v                     | £147.00 |
| 5774 | Practitioner Otoscope 3.6v Rechargeable        | £232.50 |
| 6805 | Duracell AA Size Batteries, pk/4               | £2.90   |
| 5797 | Lithium Mini Charger                           | £97.00  |
| 5813 | Bulb for 2.8v F.O./Practitioner Otoscope, pk/2 | £50.00  |
| 5814 | Bulb for 3.6v F.O./Practitioner Otoscope, pk/2 | £50.00  |
| 5791 | Disposable specula, pk/500 – specify size      | £49.00  |

## Keeler 2.8v/3.6v Fibre Optic Otoscope

Choice of 2.8v dry cell or 3.6v rechargeable lithium battery



- Genuine Halogen illumination transmitted via optic fibres for homogenous illumination leading to clear and accurate diagnosis for aural images
- Large head design for unobstructed view
- Closed head system ideally suited to pneumatic testing for tympanic mobility
- Large slide aside 2.5x magnifying lens for clear view of aural structures and wax and foreign body removal procedures under magnification.
- Supplied in a hard case with 4 x reusable specula spare bulb, pneumatic testing adaptor
- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

|      |  |         |
|------|--|---------|
| 5775 | Fibre Optic Otoscope 2.8v                      | £173.50 |
| 5776 | Fibre Optic Otoscope 3.6v Rechargeable         | £258.50 |
| 6805 | Duracell AA Size Batteries, pk/4               | £2.90   |
| 5797 | Lithium Mini Charger                           | £97.00  |
| 5813 | Bulb for 2.8v F.O./Practitioner Otoscope, pk/2 | £50.00  |
| 5814 | Bulb for 3.6v F.O./Practitioner Otoscope, pk/2 | £50.00  |
| 5791 | Disposable specula, pk/500 – specify size      | £49.00  |



# ophthalmoscopes & otoscopes

## Keeler 2.8v Jazz Diagnostic Set with Fibre Optic Otoscope



**Keeler**



**NEW**

### Ophthalmoscope Lenses



### Jazz Ophthalmoscope Head

- +/- 20D range of powers
- Renowned high performance Keeler optics with aspherical lens
- Total of 6 diagnostic beams

### Jazz Fibre Optic Otoscope Head

- 3x magnifying lens
- Brilliant Xenon illumination
- Supplied with 1 x Jazz Ophthalmoscope with handle, 1 x Jazz Otoscope with handle, 3 interchangeable colour rings in a soft case
- Requires 4 x AA batteries - not supplied
- Disposable specula available in 2.5, 3.5, 4.5 and 5.5mm

|      |   |         |
|------|---|---------|
| 9958 | Jazz Diagnostic Set 2.8v                  | £195.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 5806 | Bulb for Jazz Ophthalmoscope, pk/2        | £50.00  |
| 5811 | Bulb for Jazz Otoscope, pk/2              | £50.00  |
| 5791 | Disposable specula, pk/500 - specify size | £49.00  |

## Keeler 2.8v/3.6v Practitioner Set with Fibre Optic Otoscope



### Ophthalmoscope Lenses



### Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

#### Practitioner Ophthalmoscope Head:

- The positive maywheel lens control allows you to select a range of lens powers from +40D to -25D - increased lens range this allows detailed examination of the ocular structures and ease of use for spectacle wearers
- The Practitioner includes 6 beams for a full diagnostic capability and Blue filter, useful for assessing corneal scratches in children
- A Red-free filter for enhanced vessel examination can be used with both diagnostic beams, particularly useful for diabetic assessment and monitoring
- Soft integral brow rest allows for easy positioning and protects your lenses if you are a spectacle wearer

#### Practitioner Fibre Optic Otoscope Head:

- Genuine Halogen illumination via optic fibres for homogenous illumination leading to clear and accurate diagnosis for aural images
- Large head design for unobstructed view
- Closed head system suited to pneumatic testing for tympanic mobility
- Large slide aside 2.5x magnifying lens for clear view of aural structures and wax and foreign body removal procedures under magnification

The set is supplied in a hard case with 1 x slimline handle, 1 x replacement ophthalmoscope bulb, 1 x replacement otoscope bulb, 1 x Pneumatic testing adaptor, 5 x otoscope reusable specula and 4 x otoscope disposable specula

- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

|      |  |         |
|------|--|---------|
| 5782 | Practitioner F.O. Set, 2.8v                                | £309.00 |
| 5783 | Practitioner F.O. Set, 3.6v Rechargeable - without charger | £393.00 |
| 6805 | Duracell AA Size Batteries, pk/4                           | £2.90   |
| 5797 | Lithium Mini Charger                                       | £97.00  |
| 5807 | Bulb for 2.8v Xenon Ophthalmoscope, pk/2                   | £50.00  |
| 5813 | Bulb for 2.8v F.O./Practitioner Otoscope, pk/2             | £50.00  |
| 5808 | Bulb for 3.6v Xenon Ophthalmoscope, pk/2                   | £50.00  |
| 5814 | Bulb for 3.6v F.O./Practitioner Otoscope F.O., pk/2        | £50.00  |
| 5791 | Disposable specula, pk/500 - specify size                  | £49.00  |

# ophthalmoscopes & otoscopes

## Keeler 2.8v/3.6v Standard Set with Fibre Optic Otoscope



Keeler

*Now with Fibre Optic Otoscope*

Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

### Standard Ophthalmoscope Head:

- The positive maywheel lens control allowing you to click through lens powers from +40D to -25D allows detailed examination of the ocular structures and provides ease of use for spectacle wearers
- Two diagnostic beams - a wide angle beam for general examination and a macular beam for more detailed examination of the central retina
- Red free filter for enhanced vessel examination can be used with both diagnostic beams, and particularly useful for diabetic assessment and monitoring
- The integral soft brow rest allows for easy positioning and protects your lenses if you are a spectacle wearer

### Standard Fibre Optic Otoscope Head:

- Closed head design particularly suited to pneumatic testing
- 2.5x magnifying lens removable to allow the introduction of small instruments
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm



The set is supplied in a hard case with 1 x slimline handle, 1 x replacement ophthalmoscope bulb, 1 x replacement otoscope bulb, 1 x Pneumatic testing adaptor, 1 x bulb extractor for screw-in bulbs, 5 x otoscope reusable specula

- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied

|      |  |         |
|------|--|---------|
| 5792 | Standard Set, Dry Cell 2.8v                      | £217.00 |
| 5779 | Standard Set, 3.6v Rechargeable, without charger | £319.50 |
| 6805 | Duracell AA size Battery, pk/4                   | £2.90   |
| 5797 | Lithium Mini Charger                             | £97.00  |
| 5807 | 2.8v Xenon Ophthalmoscope Bulb, pk/2             | £50.00  |
| 5813 | Bulb for 2.8v F.O. Otoscope, pk/2                | £50.00  |
| 5808 | 3.6v Xenon Ophthalmoscope Bulb, pk/2             | £50.00  |
| 5812 | 3.6v Standard Otoscope Bulb, pk/2                | £50.00  |
| 5791 | Disposable specula, pk/500 - specify size        | £49.00  |

### Ophthalmoscope Lenses



## Keeler 2.8v/3.6v Practitioner Set with Fibre Optic Otoscope

Choice of 2.8v dry cell or 3.6v rechargeable lithium battery

### Practitioner Ophthalmoscope Head:

- The positive maywheel lens control allows you to select a range of lens powers from +40D to -25D. This allows detailed examination of the ocular structures and ease of use for spectacle wearers
- The Practitioner includes 6 beams for a full diagnostic capability and a Blue filter, useful for assessing corneal scratches in children
- Red free filter for enhanced vessel examination can be used with both diagnostic beams, is particularly useful for diabetic assessment and monitoring

### Practitioner Fibre Optic Otoscope Head:

- **Unobstructed open style head for easy** introduction of small instruments for wax or foreign body removal
- Large slide aside 2.5x lens for clear view of aural structures and procedures under magnification
- Facility to undertake pneumatic testing for tympanic mobility
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm



The set is supplied in a hard case with 1 x slimline handle, 1 x replacement ophthalmoscope bulb, 1 x replacement otoscope bulb, 1 x Airtight chamber/Pneumatic testing adaptor, 4 x otoscope disposable specula and 5 x otoscope reusable specula

- Dry cell version requires 2 x AA size batteries - not supplied
- Lithium rechargeable version requires charger - not supplied

|      |  |         |
|------|--|---------|
| 5780 | Practitioner Set, 2.8v                               | £277.00 |
| 5781 | Practitioner Set, 3.6v Rechargeable, without charger | £361.00 |
| 6805 | Duracell AA Size Batteries, pk/4                     | £2.90   |
| 5797 | Lithium Mini Charger                                 | £97.00  |
| 5807 | Bulb for 2.8v Ophthalmoscope, pk/2                   | £50.00  |
| 5813 | Bulb for 2.8v F.O./Practitioner Otoscope, pk/2       | £50.00  |
| 5808 | Bulb for 3.6v Xenon Ophthalmoscope, pk/2             | £50.00  |
| 5814 | Bulb for 3.6v F.O./Practitioner Otoscope, pk/2       | £50.00  |
| 5791 | Disposable specula, pk/500 - specify size            | £49.00  |

### Ophthalmoscope Lenses



# ophthalmoscopes & otoscopes

## Keeler Bulbs



# Keeler



- All Keeler instruments use genuine Halogen and Xenon technology, yielding a brighter, whiter illumination and a longer life
- Look at the bottom of your bulb to identify if your instrument is 2.8v dry cell battery or 3.6v rechargeable battery

### Ophthalmoscope Bulbs

|      |   |          |        |
|------|---|----------|--------|
| 5806 | Pocket Ophthalmoscope Bulb 2.8v Halogen, pk/2           | <b>1</b> | £50.00 |
| 5807 | Standard, Practitioner Ophthalmoscope 2.8v Bulbs, pk/2  | <b>2</b> | £50.00 |
| 5808 | Standard, Practitioner Ophthalmoscope 3.6v Bulbs, pk/2  |          | £50.00 |
| 5809 | Vista & Vista 20 Ophthalmoscope Bulb 2.8v Halogen, pk/2 | <b>3</b> | £50.00 |
| 5810 | Vista & Vista 20 Ophthalmoscope Bulb 3.6v Halogen, pk/2 |          | £50.00 |

### Otoscope Bulbs - screw in

|      |   |          |        |
|------|---|----------|--------|
| 5811 | Standard, Pocket Deluxe, Vetscope Otoscope & Tongue Depressor Bulb - 2.8v Halogen, pk/2 | <b>4</b> | £50.00 |
| 5812 | Standard, Deluxe, Vetscope Otoscope and Tongue Depressor Bulb - 3.6v Halogen, pk/2      |          | £50.00 |

### Otoscope Bulbs - push in

|      |   |          |        |
|------|---|----------|--------|
| 5813 | Practitioner, Vista, Fibre Optic Otoscope & Finhoff Transilluminator Bulb, 2.8v Halogen, pk/2   | <b>5</b> | £50.00 |
| 5814 | Practitioner, Vista, Fibre Optic Otoscope and Finhoff Transilluminator Bulb, 3.6v Halogen, pk/2 |          | £50.00 |

## Keeler Otoscope Disposable Specula


- Suitable for Practitioner, Standard, Pocket, Vista, Fibre Optic



|         |                                   |        |
|---------|-----------------------------------|--------|
| 5791-25 | Disposable specula, 2.5mm, pk/500 | £49.00 |
| 5791-35 | Disposable specula, 3.5mm, pk/500 | £49.00 |
| 5791-45 | Disposable specula, 4.5mm, pk/500 | £49.00 |
| 5791-55 | Disposable specula, 5.5mm, pk/500 | £49.00 |

## Keeler Otoscope Reusable Specula

- Suitable for Practitioner, Vista, Fibre Optic, Standard, Pocket



|      |                         |        |
|------|-------------------------|--------|
| 5790 | Reusable specula, set/5 | £25.00 |
|------|-------------------------|--------|

## Keeler Lithium Mini Charger



- Fast charges in 2 hours
- Provides over 2 hours of continuous use
- Can use instrument when charging
- LEDs indicate when battery life is low, charging or fully charged

|      |   |        |
|------|---|--------|
| 5797 | Lithium Mini Charger for Slim Line Handle | £97.00 |
|------|---|--------|

## Keeler Lithium Dual Charger



- Recharges two Slim Line handles
- Fast charges in 2 hours of continuous use
- Provides over 2 hours use time
- LEDs indicate when battery life is low, charging or fully charged

|      |                                      |         |
|------|--------------------------------------|---------|
| 5798 | Lithium Charger for Slim Line handle | £212.00 |
|------|--------------------------------------|---------|

## Keeler General Accessories



|      |  |          |        |
|------|--|----------|--------|
| 5822 | Bulb Extractor for Both Old & New Type Standard & Pocket             | <b>6</b> | £3.90  |
| 5796 | Airtight Chamber/Magnifier for Both Old & New Type Standard & Pocket | <b>7</b> | £16.00 |
| 5821 | Otoscope Viewing Lens, for Practitioners & Vista                     |          | £15.50 |

## Keeler 3.75v Lithium Battery



- Lithium battery for use with all Keeler slim line instrument handles
- When fully charged provides over 2 hours use time

|      |   |        |
|------|---|--------|
| 5805 | Rechargeable Lithium Battery for Slim Line Handle | £75.00 |
| 5800 | Red End Cap for Slim Line Handle                  | £24.00 |

# ophthalmoscopes & otoscopes

## Keeler 3.6v Rechargeable Desk Set



**Keeler**

- 3.6v Illumination
- Can be desk or wall mounted to suit your needs
- Lithium Ion technology provides over 2 hour of continuous use and 2 hour rapid charge
- LEDs indicate when battery life is low, charging or fully charged
- Complete with two handles and charger
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

### Choice of 3 combinations

- Set A Standard Ophthalmoscope, Standard Otoscope
- Set B Practitioner Ophthalmoscope, Practitioner Otoscope
- Set C Practitioner Ophthalmoscope, Practitioner Fibre Optic Otoscope

|      |   |          |         |
|------|---|----------|---------|
| 5784 | Rechargeable Set A Standard Set                               | <b>1</b> | £635.00 |
| 5785 | Rechargeable Set B Practitioner Set                           |          | £669.00 |
| 5786 | Rechargeable Set C Practitioner with Fibre Optic Otoscope Set |          | £756.00 |
| 5799 | Charger Wall Mounting Kit                                     |          | £24.00  |
| 5791 | Disposable specula, pk/500 - specify size                     |          | £49.00  |



## Keeler 3.6v Cord Handle Wall Set

- 3.6v Illumination
- Keeler's anti-theft designed ophthalmoscope and otoscopes have a lockable head and handle fitting
- Strong coiled tangle-free cable permits instrument use up to 3 metres from transformer position
- Automatic illumination when the instruments are removed from their cradle. With LED confirmation
- Mains only
- Disposable specula available in 2.5mm, 3.5mm, 4.5mm and 5.5mm

### Choice of 3 combinations

- Set A Standard Ophthalmoscope, Standard Otoscope
- Set B Practitioner Ophthalmoscope, Practitioner Otoscope
- Set C Practitioner Ophthalmoscope, Practitioner Fibre Optic Otoscope

- Optional wall mounted Dispos-A-Spec, specula dispenser supplied with 25 x 2.5 mm, 3.5 mm, 4.5 mm and 5.5 mm disposable specula

|      |   |          |         |
|------|---|----------|---------|
| 5787 | Cord Handle Set A Standard Set                        |          | £646.00 |
| 5788 | Cord Handle Set B Practitioner Set                    |          | £673.00 |
| 5789 | Cord Handle Set C Practitioner with F.O. Otoscope set | <b>2</b> | £765.00 |
| 5820 | Dispos-A-Spec   | <b>3</b> | £61.00  |
| 5791 | Disposable specula, pk/500 - specify size             |          | £49.00  |



# ophthalmoscopes & otoscopes

## Heine 2.5v K180 Fibre Optic Set



Heine K180 Diagnostic Set comprises:

### K180 Otoscope Head

- 2.5v fibre optic illumination
- High quality acrylic viewing window with 3x magnification
- High output Xenon bulb
- Swivelling viewing window
- Insufflation port
- Supplied with 4 x reusable specula
- Disposable specula available in 2.5mm and 4.0mm

### K180 Direct Ophthalmoscope Head:

- Soft orbital rest, ergonomical shape
- 5 apertures; slit, red-free, target fixation, large and small spot
- Aspherical optics give brilliant image with reduced corneal reflexes
- Dust-resistant housing
- 28 lenses from - 35 D to + 40 D

- Supplied with metal handle in case
- Requires 2 x C size batteries – not supplied

### Ophthalmoscope Lenses



|      |   |         |
|------|---|---------|
| 6967 | K180 Fibre Optic Set                      | £269.00 |
| 6806 | Duracell C Size Battery, pk/2             | £2.90   |
| 7540 | Bulb 037 for 2.5v Otoscope                | £16.00  |
| 7539 | Bulb 184 for 2.5v Ophthalmoscope          | £16.00  |
| 7118 | Disposable specula, pk/100 - specify size | £10.00  |

## Heine 2.5v K100 Otoscope



- 2.5v direct illumination
- High quality acrylic viewing window with 3x magnification
- High output Xenon bulb
- Swivelling viewing window
- Insufflation port
- Metal handle
- Supplied with 3 x reusable specula in soft canvas case
- Requires 2 x C size batteries – not supplied
- Disposable specula available in 2.5mm and 4.0mm

|      |   |         |
|------|---|---------|
| 6978 | K100 Direct Illumination Otoscope         | £128.00 |
| 6806 | Duracell C Size Battery, pk/2             | £2.90   |
| 7540 | Bulb 037 for 2.5v D.I. Otoscope Bulb      | £16.00  |
| 7583 | Disposable specula, pk/100 - specify size | £10.00  |

# ophthalmoscopes & otoscopes

## Heine Fibre Optic Modular Wall System



### K Series Wall Unit

- Modular Wall unit supplied as standard with 2 instrument bays
- Heine 3.5v K180 Ophthalmoscope, lens range +20D to -35D
- Heine 3.5v K100 Diagnostic Otoscope or K180 Fibre Optic Otoscope
- 3m spiral cord handles with magnetic automatic on/off and ceramic rheostats
- Instrument locks as standard to prevent unauthorised removal
- Disposable and reusable tips
- Hygienic easy to clean wall mounted unit, smooth rounded corners, no hidden recesses

#### Optional Modular Equipment:

- 3rd instrument bay, with recommended additional instrument 3.5v tongue depressor with blade ejector and 100 disposable blades
- All Tip Specula dispenser
- XXL large dial aneroid with choice of 1, 2 or 3 cuff holders

#### Ophthalmoscope Lenses



### Module 1



Module 1, the principal element, choice between K or Beta Series instruments.

6968 K Series Wall Unit

£577.00

### Beta Series Wall Unit

- Modular Wall unit supplied as standard with 2 instrument bays
- Heine 3.5v B200 Ophthalmoscope, lens range +20D to -35D
- Heine 3.5v B100 Diagnostic Otoscope or B200 Fibre Optic Otoscope
- 3m spiral cord handles with magnetic automatic on/off and ceramic rheostats
- Instrument locks as standard to prevent unauthorised removal
- Disposable and reusable tips
- Hygienic easy to clean wall mounted unit, smooth rounded corners no hidden recesses

#### Optional Modular Equipment:

- 3rd instrument bay, with recommended additional instrument 3.5v tongue depressor with blade ejector and 100 disposable blades
- All Tip Specula dispenser
- XXL large dial aneroid with choice of 1, 2 or 3 cuff holders

#### Ophthalmoscope Lenses



6969 Beta Series Wall Unit

£665.00

## Optional Equipment Modules 2, 3 and 4 for Heine Modular Wall System above

Module 2 Instrument Bay, code 6970. Illustration shows the 3.5v tongue depressor with blade ejector code 6971 fitted to the handle, supplied with 100 disposable blades

|      |   |   |         |
|------|---|---|---------|
| 6970 | 3rd Instrument Bay, excludes instrument head                        | 2 | £159.00 |
| 6971 | 3.5v tongue depressor with blade ejector head for disposable blades | 3 | £66.35  |
| 6972 | 3.5v tongue depressor with blade ejector head for wooden blades     |   | £66.35  |
| 6973 | Disposable Blades, pk/100   |   | £8.35   |

Module 3 Heine Allspec Tip Dispenser, the convenient and hygienic way to store and dispense otoscope specula. Capacity 200 specula; 100 x 2.5mm and 100 x 4.0mm size


|      |                             |   |        |
|------|-----------------------------|---|--------|
| 6975 | Heine Allspec Tip Dispenser | 4 | £24.00 |
|------|-----------------------------|---|--------|

#### Heine Otoscope specula

| SAVE MORE |                                  | 1      | 5+    |
|-----------|----------------------------------|--------|-------|
| 7118-25   | Disposable specula 2.5mm, pk/100 | £10.00 | £8.00 |
| 7118-40   | Disposable specula 4.0mm, pk/100 | £10.00 | £8.00 |

Module 4 Gamma XXL large dial aneroid sphygmomanometer with adult latex free cuff

|      |  |   |        |
|------|--|---|--------|
| 6974 | Gamma XXL Sphyg, with Adult Cuff       | 5 | £99.00 |
| 7124 | Small Child Cuff 8-13cm for Gamma XXL  |   | £20.30 |
| 7125 | Child Cuff 13-20cm for Gamma XXL       |   | £26.90 |
| 7126 | Small Adult Cuff 20-29cm for Gamma XXL |   | £26.90 |
| 7127 | Adult Cuff 29-41cm for Gamma XXL       |   | £27.40 |
| 7128 | Thigh Cuff up to 41cm for Gamma XXL    |   | £35.60 |



Module 2 Instrument Bay

Module 4 Gamma XXL



# ophthalmoscopes & otoscopes

## Heine 2.5v M3000 Otoscopes



- Modern 2.5v pocket otoscope, unique compact design
- Enhanced XHL Xenon Halogen Technology with 100% more light
- Viewing window with 3x magnification
- Swivelling viewing window built into instrument
- Pocket clip with integrated on/off switch. Switches off automatically when returned to the pocket
- High quality shockproof handle, chrome-finish upper section
- Supplied in case with 4 x reusable specula and 5 x 2.5mm disposable specula & 5 x 4mm disposable specula
- Requires 2 x AA batteries – not supplied
- Disposable specula available in 2.5mm and 4.0mm

|      |   |        |
|------|---|--------|
| 6976 | M3000 Otoscope                            | £87.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90  |
| 7117 | Bulb 110 for 2.5v M3000 Otoscope          | £16.00 |
| 7118 | Disposable specula, pk/100 - specify size | £10.00 |

## Heine 2.5v M3000 Ophthalmoscope



- Modern 2.5v pocket ophthalmoscope, unique compact design
- Outstanding optics provide crisp retinal detail
- -20 to +20 with 18 lenses, -/+1 to -4/+4 in single steps
- 5 apertures for a full diagnostic capability, large spot, small spot, fixation star for eccentric fixation, hemispot for small pupils, red-free filter
- Pocket clip with integrated on/off switch. Switches off automatically when returned to the pocket
- High quality shockproof handle, chrome-finish upper section
- Supplied in hard case
- Requires 2 x AA batteries – not supplied

### Ophthalmoscope Lenses



|      |  |         |
|------|--|---------|
| 6979 | M3000 Ophthalmoscope                   | £139.00 |
| 6805 | Duracell AA Size Batteries, pk/ 4      | £2.90   |
| 7116 | Bulb 106 for 2.5v M3000 Ophthalmoscope | £16.00  |

## Heine 2.5v M3000 Set



- Set comprises of:  
M3000 Otoscope with 2.5v handle  
M3000 Ophthalmoscope with 2.5v handle
- Supplied in case with 4 x reusable specula and 5 x 2.5mm disposable specula & 5 x 4mm disposable specula
- Each handle requires 2 x AA batteries – not supplied
- Disposable specula available in 2.5mm and 4.0mm

### Ophthalmoscope Lenses



|      |   |         |
|------|---|---------|
| 6965 | M3000 2 Handle Set                        | £211.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 7117 | Bulb 110 for 2.5v M3000 Otoscope          | £16.00  |
| 7116 | Bulb 106 for 2.5v M3000 Ophthalmoscope    | £16.00  |
| 7118 | Disposable specula, pk/100 - specify size | £10.00  |

## Heine 2.5v M3000 Multi Set



- Set comprises of:  
M3000 2.5v handle  
M3000 Direct Illumination Otoscope head  
M3000 Ophthalmoscope head  
M3000 Laryngeal mirror, straight stem  
M3000 Throat illuminator and tongue depressor
- Supplied in hard case with 4 x reusable specula
- Requires 2 x AA batteries – not supplied
- Disposable specula available in 2.5mm and 4.0mm

### Ophthalmoscope Lenses



|      |   |         |
|------|---|---------|
| 6966 | M3000 Multi Set                           | £239.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 7117 | Bulb 110 for 2.5v M3000 Otoscope          | £16.00  |
| 7116 | Bulb 106 for 2.5v M3000 Ophthalmoscope    | £16.00  |
| 7118 | Disposable specula, pk/100 - specify size | £10.00  |

# ophthalmoscopes & otoscopes

## Riester pen-scope® 2.7v Pocket Otoscope & Ophthalmoscope

**Riester**  
The familiar way

- Practical metal clip on the handle for attaching to the doctor's coat
- Otoscope with direct illumination, 2.5x magnifying lens, removal for inserting external instruments
- Ophthalmoscope with four different apertures and correction lenses from 0 to +20 and 0 to -20 diopters, easy to operate focusing wheel and integrated eyeglass protector
- Available in Black, Saffron Yellow or Green, specify colour when ordering
- Requires 2 x AA size batteries - not supplied
- Disposable or reusable specula available in 2.5mm and 4.0mm

|      |  |         |
|------|--|---------|
| 5875 | pen-scope® Otoscope – specify colour       | £57.00  |
| 5879 | pen-scope® Ophthalmoscope – specify colour | £91.00  |
| 5883 | pen-scope® Otoscope & Ophthalmoscope Set   | £136.00 |
| 6805 | Duracell AA Size Batteries, pk/4           | £2.90   |
| 7579 | Disposable specula, pk/500 - specify size  | £51.30  |
| 5896 | Reuseable specula, pk/10 - specify size    | £17.00  |
| 5902 | Bulb for pen-scope® Otoscope               | £22.80  |
| 7565 | Bulb for pen-scope® Ophthalmoscope         | £22.80  |



Ophthalmoscope Lenses



## Riester ri-mini® 2.5v Pocket Otoscope & Ophthalmoscopes

### Fibre Optic Otoscope

- Magnifying lens with 2.5x magnification, removable

### Ophthalmoscope

- Five different apertures and correction lenses from 0 to +20 and 0 to -20 diopters, easy to operate focusing wheel and integrated eyeglass protector

### Set

- Metal clip on the handle
- Supplied in protective hard case as individual instruments or as a two handle set
- Available in Black, Saffron Yellow or Green, specify colour when ordering
- Requires 2 x AA size batteries - not supplied
- Disposable or reusable specula available in 2.5 and 4.0mm

|      |   |         |
|------|---|---------|
| 5858 | ri-mini® Otoscope – specify colour                      | £79.50  |
| 5865 | ri-mini® Ophthalmoscope – specify colour                | £114.00 |
| 5871 | ri-mini® Otoscope & Ophthalmoscope Set – specify colour | £171.00 |
| 6805 | Duracell AA Size Batteries, pk/4                        | £2.90   |
| 7579 | Disposable specula, pk/500 - specify size               | £51.30  |
| 5896 | Reuseable specula, pk/10 - specify size                 | £17.00  |
| 5901 | Bulb for ri-mini® Otoscope                              | £22.80  |
| 5903 | Bulb for ri-mini® Ophthalmoscope                        | £22.80  |



Ophthalmoscope Lenses



## Riester e-scope® LED Fibre Optic Otoscope & Ophthalmoscope

### e-scope® LED Fibre Optic Otoscope

- Triple magnification for optimum view
- Sturdy ear speculum receptacle in hardwearing hygienic metal
- Suitable for performance of pneumatic tests - supplied without connector and insufflators
- Disposable specula available in 2.5mm and 4.0mm

### e-scope® LED Ophthalmoscope

- 3.7v LED illumination for higher-contrast and more efficient diagnosis
- Dioptre disc with 18 corrective lenses - +/- 1, 2, 3, 4, 6, 8, 10, 15, 20 dioptres
- Total of 6 diagnostic beams
- High performance optics with aspherical condenser lens
- Spectacle protection
- Parallel optical path
- Dust protected

- Supplied as individual instruments or as a two handle set in a case
- Requires 2 x AA size batteries - not supplied

|      |   |         |
|------|---|---------|
| 9961 | e-scope® LED Fibre Optic Otoscope         | £85.00  |
| 9949 | e-scope® LED Ophthalmoscope               | £150.00 |
| 9950 | e-scope® LED Set                          | £225.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 7579 | Disposable specula, pk/500 - specify size | £51.30  |



Ophthalmoscope Lenses



# ophthalmoscopes & otoscopes

## Riester ri-scope® 2.5v L1 Direct Illumination and L2 Fibre Optic Otosopes



- L1 – 2.5v halogen direct illumination otoscope
- L2 – 2.5v fibre optic otoscope
- Bi-directional swivelling optical glass with 3x magnification
- Suitable for pneumatic tests – supplied without ball
- C-Cell battery handle
- Requires 2 x AA size batteries – not supplied
- Disposable or reusable specula available in 2.5mm and 4.0mm



## Riester ri-scope® L3 Fibre Optic Otosopes

**Riester**  
The familiar way



- Available as 2.5v Halogen or 3.5v LED illumination
- Fibre optics for optimal beam and light transmission
- Bi-directional swivelling optical glass with 3x magnification
- Operation lens with 4x magnification and glass mini optic ease the insertion of instruments
- Easy to operate SpecEjec specula ejection device
- Fitting for ri-scope® specula and enlarged conduit for instruments
- Suitable for pneumatic tests – supplied without ball
- Disposable or reusable specula available in 2.0mm and 4.0mm



|      |   |         |
|------|---|---------|
| 5846 | L1 Direct Illumination 2.5v Otscope       | £119.00 |
| 5847 | L2 Fibre Optic 2.5v Otscope               | £148.00 |
| 6805 | Duracell AA Size Batteries, pk/4          | £2.90   |
| 7579 | Disposable specula, pk/500 – specify size | £51.30  |
| 5896 | Reusable specula, pk/10 – specify size    | £17.00  |
| 5898 | Bulb for L1 Otscope                       | £22.80  |
| 5901 | Bulb for L2/L3 Otscope                    | £22.80  |

|      |   |         |
|------|---|---------|
| 9439 | L3 2.5v HL Otscope                        | £176.70 |
| 9440 | L3 3.5v LED Otscope                       | £233.70 |
| 9446 | Disposable specula, pk/500 – specify size | £51.30  |
| 9444 | Reusable specula, pk/10 – specify size    | £17.00  |
| 5901 | Bulb for 2.5v L3 Otscope                  | £22.80  |
| 9452 | Bulb for 3.5v LED L3 Otscope              | £22.80  |
| 9443 | Ri-charger® for 3.5v handles              | £171.00 |

## Riester ri-scope® 2.5v L1 and L2 Ophthalmoscopes



- L1 – with 7 different apertures
- L2 – with 9 different apertures
- 2.5v halogen illumination
- Diopter wheel with 29 corrective lenses:  
Plus 1-10, 12, 15, 20, 40  
Minus 1-10, 15, 20, 25, 30, 35
- Easy to operate aperture hand-wheel
- Optimised high performance optics with aspherical condenser lens
- Parallel beam path
- Dust proof
- C-Cell battery handle
- Requires 2 x AA size batteries – not supplied



Ophthalmoscope Lenses



|      |                                  |         |
|------|----------------------------------|---------|
| 5848 | L1 2.5v Ophthalmoscope           | £171.00 |
| 5849 | L2 2.5v Ophthalmoscope           | £193.00 |
| 6805 | Duracell AA Size Batteries, pk/4 | £2.90   |
| 5903 | Bulb for L1/L2/L3 Ophthalmoscope | £22.80  |

## Riester ri-scope® L3 Ophthalmoscopes



- Available as 2.5v Halogen or 3.5v Xenon versions
- Corrective lenses for 89 diopter values – Plus 1-45 in single steps and Minus 1-44 in single steps
- Easy to operate aperture hand wheel with 10 different apertures
- Focusing device with zero retention force for the quick correction of ametropias
- Dust proof



Ophthalmoscope Lenses



|      |                                 |         |
|------|---------------------------------|---------|
| 9441 | L3 2.5v HL Ophthalmoscope       | £205.20 |
| 9442 | L3 3.5v XL Ophthalmoscope       | £290.70 |
| 5903 | Bulb for 2.5v L3 Ophthalmoscope | £22.80  |
| 9453 | Bulb for 3.5v L3 Ophthalmoscope | £22.80  |
| 9443 | Ri-charger® for 3.5v handles    | £171.00 |

# ophthalmoscopes & otoscopes

## Riester ri-former® 3.5v Modular Wall Diagnostic Station

**Riester**  
The familiar way



- Premium diagnostic station with 3.5v L3 Ophthalmoscope, Otoscope, Big Ben Sphygmomanometer, Digital clock and specula dispenser
- Anti-theft to ensure your instruments are where you left them
- Non-slip, ergonomically designed handles with rheostat, the handle enables continuously adjustable brightness as well as manual switching on and off
- Spiral cord stretches up to 3m
- Easy cleaning and disinfection of the instruments
- Optional ri-thermo® infra red thermometer with anti-theft security
- Disposable or reusable specula available in 2.0mm and 4.0mm

|      |   |         |
|------|---|---------|
| 9449 | L3 3.5v Wall Diagnostic Station                 | £875.00 |
| 9450 | Ri-thermo® infra red thermometer                | £205.20 |
| 9451 | Probe covers for ri-thermo® thermometer, pk/100 | £22.80  |
| 9446 | Disposable specula, pk/500 – specify size       | £51.30  |
| 9444 | Reusable specula, pk/10 – specify size          | £17.00  |

**2 Year Warranty**

## Riester 2.5v L1 or L2 Wall Diagnostic Stations



Optional wall mounted Big Ben Sphyg with adult cuff

- Available in two versions; with the L1 otoscope and ophthalmoscope heads or the L2 otoscope and ophthalmoscope heads and digital clock
- Anti-theft to ensure your instruments are where you left them
- Optional ri-thermo® infra red thermometer with anti-theft security
- Optional ear specula dispenser
- Optional Big Ben aneroid wall mounting sphygmomanometer and ear specula dispenser available to allow you the flexibility to build a comprehensive wall diagnostic station
- Disposable or reusable specula available in 2.5mm and 4.0mm

|      |  |          |         |
|------|--|----------|---------|
| 5855 | L1 Wall Unit                                   | <b>1</b> | £568.00 |
| 5854 | L2 Wall Unit                                   |          | £627.00 |
| 5856 | Big Ben Sphygmomanometer, wall mounted         | <b>2</b> | £79.00  |
| 5857 | Riester Ear specula Dispenser, wall mounted    | <b>3</b> | £68.00  |
| 9451 | Ri-thermo infra red thermometer                |          | £205.20 |
| 9450 | Probe covers for ri-thermo thermometer, pk/100 |          | £22.80  |
| 7579 | Disposable specula, pk/500 - specify size      |          | £51.30  |
| 5896 | Reusable specula, pk/10 - specify size         |          | £17.00  |

### Replacement Bulbs

|      |                               |        |
|------|-------------------------------|--------|
| 5903 | Bulb for L1/L2 Ophthalmoscope | £22.80 |
| 5898 | Bulb for L1 Otoscope          | £22.80 |
| 5901 | Bulb for L2 Otoscope          | £22.80 |



Optional Ear specula Dispenser complete with starter pack of specula



**2 Year Warranty**



# ophthalmoscopes & otoscopes

## Riester 2.5v L1/L2 Sets



**Riester**  
The familiar way

- 2.5v halogen illumination
- Available in two versions; with the L1 or L2 otoscope and ophthalmoscope heads, single handle
- Supplied case with starter pack of otoscope specula, spare otoscope and ophthalmoscope bulbs
- Disposable or Reusable specula available in 2.5mm and 4.0mm

|      |                                  |         |
|------|----------------------------------|---------|
| 5850 | ri-scope® L1 Set                 | £256.50 |
| 5851 | ri-scope® L2 F.O. Set            | £302.00 |
| 6805 | Duracell AA Size Batteries, pk/4 | £2.90   |

### Specula

|      |   |        |
|------|---|--------|
| 7579 | Disposable specula, pk/500 - specify size | £51.30 |
| 5896 | Reusable specula, pk/10 - specify size    | £17.00 |

### Replacement Bulbs

|      |                               |        |
|------|-------------------------------|--------|
| 5903 | Bulb for L1/L2 Ophthalmoscope | £22.80 |
| 5898 | Bulb for L1 Otoscope          | £22.80 |
| 5901 | Bulb for L2 Otoscope          | £22.80 |

## Riester 3.5v Rechargeable Desk Units



**Riester**  
The familiar way

- 3.5v LED illumination in L2 otoscopes provides brighter longer lasting bulb
- 3.5v Xenon illumination
- Available in two versions; with the L1 or L2 otoscope and ophthalmoscope heads
- Indicator lights show charging yellow light, ready for use green light
- Can be wall mounted
- Centre compartment to place otoscope specula
- Supplied with 3.5v Lithium rechargeable batteries
- Disposable or Reusable specula available in 2.5mm and 4.0mm

|      |                                |         |
|------|--------------------------------|---------|
| 5852 | L1 Rechargeable Desk Unit      | £507.00 |
| 5853 | L2 F.O. Rechargeable Desk Unit | £598.50 |

### Specula

|      |   |        |
|------|---|--------|
| 7579 | Disposable specula, pk/500 - specify size | £51.30 |
| 5896 | Reusable specula, pk/10 - specify size    | £17.00 |

### Replacement Bulbs

|      |                               |        |
|------|-------------------------------|--------|
| 5903 | Bulb for L1/L2 Ophthalmoscope | £22.80 |
| 5898 | Bulb for L1 Otoscope          | £22.80 |
| 5901 | Bulb for L2 Otoscope          | £22.80 |

## Dinolite Digital USB Video Otoscope



- High resolution Digital Video Otoscope
- Simply plugs into any available USB port for an instant high resolution image of your patient's ear canal
- Built-in light source and camera give live views with option to take snapshots and video
- Software included
- Software enables annotation of images and direct e-mail of images or video files
- 5mm and 3mm disposable speculae available
- Supplied with cable, pouch, software CD, pack of wipes, 100 x 5mm disposable specula, desktop holder and full instructions

|      |   |         |
|------|---|---------|
| 9418 | Digital USB Video Otoscope                | £325.00 |
| 9419 | Disposable Specula, pk/100 - specify size | £12.95  |