

EasyOn-PC Spirometer



Connects directly into a PC or laptop for real time spirometry diagnostics

- Bio calibration check feature
- Inspiratory and expiratory real-time curve
- Test selection with or without tidal breathing
- Pre/post measurements
- Quick test
- Child incentive
- Trending
- Actualised Fletcher curve on report
- Country specific predicted values selectable
- Interpretation mode selectable
- Lung age
- Auto QC
- Ethnic correction factors
- Multilingual Modes, parameters and standards:
- Auto QC user selection for best trial from 8
- Tidal breathing selectable
- FEV1, FEV6, FEV1/FEV6, FVC, FVL
- FEV1/FVC, FEV1/VC, FEV1/VCin, PEF, Pre/Post, % var, MEF (FEF)
- ERS, ATS
- No sensor in contact with sample
- No influence of humidity, barometric pressure, contamination
- No calibration required, calibration can be verified with 3 Litre Syringe
- No cleaning, no maintenance
- Extremely high accuracy for low flows
- Weight: 155g
- Cable length: 1.8m
- Windows 2000 SP4, Windows Server 2003, Windows XP SP2, Windows Vista
- RAM: min 512MB - 1GB recommended
- Supplied with mouthpiece adaptor to allow One Way Spirometry mouthpieces to be used instead of Spirettes



Laptop Not Included

9623	EasyOn-PC Spirometer	£995.00
9624	Spirettes, pk/50	£65.00

Universal 3 Litre Calibration Syringe



- For use in verifying calibration of Spirometers
- Manufactured to fit directly onto Carefusion, Vitalograph and MIR Spirometers without any adaptor
- Supplied with silicone adaptor to allow the EasyOn-PC Spirometer

9630	Calibration Syringe	£249.00
------	---------------------	---------

Single use Spirette™ for EasyOn-PC



- Singularly wrapped for single patient use
- Bacteria filter membrane
- Suits all ages and all levels of fitness
- Smooth non-stick polyethylene
- Shaped oval for optimal lip seal
- Environmentally friendly, disposable and biodegradable

9624	Spirette™, pk/50	£65.00
------	------------------	--------

Disposable Mouthpieces



- Disposable One Way valve coated mouthpieces
- Will fit most makes of spirometers and peak flow meters
 - Prevents the patient from inhaling whilst performing a spirometry or peak flow test
 - Standard 30mm

9628	One Way Mouthpieces, pk/200	£25.00
------	-----------------------------	--------

- Disposable cardboard coated mouthpieces
- Will fit most makes of spirometers and peak flow meters
 - Standard 30mm



9629	Mouthpieces, pk/500	£23.00
------	---------------------	--------

Welch Allyn PC-Based SpiroPerfect™ Spirometer


WelchAllyn®

emis
accredited partner

systemone
ONE PATIENT ONE RECORD



Comparison Analysis



Trend Analysis

- Performs FVC, Relaxed VC, FEV1, FEV1%, PEF and pre/post testing with automatic interpretation, comparison and lung age calculation
- View real-time Flow/Volume Loop and Volume/Time graphs
- Choice of reports - on screen/printed
- Able to compare previous tests
- Technical Support Packages available
- Supplied complete with sensor, pressure tubing, calibration syringe, 4 flow transducers, CardioPerfect Workstation software, two training posters and quick reference card
- See page 86 for ECG Spiro combination

6164	Spirometer with Calibration Syringe	£1281.00
6165	Disposable Spiro Flow Transducers, pk/25	£24.00
6166	Disposable Spiro Flow Transducers, pk/100	£93.00

FREE
1 Box of 100
Transducers
worth £93.00

Demo Available

Free Welch Allyn DS-55 Durashock on Demo request until 31st March 2011

Free Welch Allyn Diagnostic set with every purchase until 31st March 2011



MicroLoop Spirometer with Gold Standard Turbine



CareFusion



systemone
DEPARTMENTAL RECORD



- Fully compliant to ATS/ERS Task Force "Standardisation of Lung Function Testing", Standardisation of Spirometry, European Respiratory Journal 2005: 26: 319-338
- High definition colour touch-screen
- Fast on-screen textual entry enabling comments and report writing
- Desktop docking station for battery charging and PC / Printer connection
- Direct connection to external Hewlett Packard printers with fully customisable printout format
- Unique contextual Help Screens (language selectable) with Help Button
- 41 inspiratory and expiratory indices measurable
- On-screen graphical child incentive devices
- Large 2000 patient test memory
- Optional Spirometry PC Software (SPCS)
- Internationally recognised on-screen test quality assurance prompts, e.g. slow start, abrupt end etc. ensuring quality test results
- 'Today's testing feature' - creating a list of patients tested or to be tested
- Numerous predicted values and languages selectable
- 'Lung Age' assessment and textual result interpretation by ATS, BTS-NICE or Enright
- Body Mass Index and Dyspnoea Score
- Bronchodilator and/or Steroid assessment
- Supplied complete with all accessories in a sturdy carrying case
- Dimensions: 12.3H x 8.1W x 2.3D cm Transducer 5H x 6W x 9D cm
- Weight: Excluding transducer Micro Loop unit 191g, Docking Station 124g

7156	MicroLoop Spirometer	£1127.50
7157	MicroLoop Spirometer with Spirometry PC Software	£1333.50

Mouthpieces

7165	Adult Disposable, pk/500	£25.00
7170	Adult Disposable One Way Valve, pk/250	£35.00

Micro and Micro Plus Spirometers



The Micro Spirometers are specifically designed for situations where low cost precision respiratory measurement is required.

- Choice of two versions;
The Micro Spirometer measures FEV1 and FVC
The MicroPlus Spirometer will make the additional measurements of PEF and FER (FEV1%) and is also fully compatible with Spirometry PC Software (SPCS)
- Supplied with hard carry case and full instructions
- Dimensions: 17H x 6W x 2.6D cm
- Weight: 190g

7151	Micro Spirometer	£200.00
7152	Micro Plus Spirometer	£303.50
7012	Spirometry PC Software (SPCS)	£395.00
7174	PC Connection Cable for use with Micro Plus	£16.50

Mouthpieces

7165	Adult Disposable, pk/500	£25.00
7170	Adult Disposable One Way Valve, pk/250	£35.00

MicroGP Spirometer



- The MicroGP is the device of choice for COPD and Asthma screening offering a high degree of accuracy with simplicity of operation
- Measuring and displaying FEV1/FVC/PEF with percent predicted as well as allowing a direct printout which includes the small airways flow indicators of F50/F25/MEF and the Flow/Volume and Volume/Time curves
- Fully compatible with Spirometry PC Software (SPCS)
- Supplied with hard carry case and full instructions
- Dimensions: 17H x 7W x 6D cm
- Weight: 190g

7153	MicroGP Spirometer	£406.50
7012	Spirometry PC Software (SPCS)	£406.50
7174	PC Connection Cable	£16.50

Mouthpieces

7165	Adult Disposable, pk/500	£25.00
7170	Adult Disposable One Way Valve, pk/250	£35.00

MicroLab Spirometer with Gold Standard Turbine



systemone
DEPARTMENT RECORD



- Fully compliant to ATS/ERS Task Force "Standardisation of Lung Function Testing", Standardisation of Spirometry, European Respiratory Journal 2005: 26: 319-338
- High definition colour touch-screen
- Simple screen icon interface using stylus or optional MicroMouse
- Fast on-screen textual entry enabling comments and report writing
- Fast, quiet printer with easy load paper mechanism
- Fully customisable printout format
- Direct connection to external Hewlett Packard printers
- 41 inspiratory and expiratory indices measurable
- On-screen graphical child incentive device
- Large 2000 patient test memory
- Optional Spirometry PC Software (SPCS)
- Internationally recognised on-screen test quality assurance prompts, e.g. slow start, abrupt end etc. ensuring quality test results
- Reversible test can be done on Asthma patients
- Unique 'today's testing feature' - creating a list of patients tested or to be tested
- Numerous predicted values and languages selectable
- 'Lung Age' assessment and textual result interpretation by ATS, BTS-NICE or Enright
- Body Mass Index and Dyspnoea Score
- Bronchodilator and/or Steroid assessment
- Supplied complete with all accessories in a sturdy carrying case
- Dimensions: 25.5H x 12W x 3.5D cm Transducer 5H x 6W x 9D cm
- Weight: Excluding any transducers 630g

7159	MicroLab Spirometer	£1333.50
7160	MicroLab Spirometer with Spirometry PC Software	£1539.50
7168	Printer paper, pk/5	£22.50

Mouthpieces

7165	Adult Disposable, pk/500	£25.00
7170	Adult Disposable One Way Valve, pk/250	£35.00

Micro Medical Accessories

7173	Nose Clip, pk/5	£3.00
7164	3 litre Calibration Syringe	£303.50



SpiroUSB



PC Not Included



CareFusion



systemone
ONEPATIENTONERECORD



- Advanced Spirometer designed for the professional that plugs directly into your PC, includes Spirometry PC Software
- Up to 41 spirometry indices are measurable, open and closed loop test techniques and the overlaying of previous test curves is possible
- Extensive search facilities and powerful trending of results together with lung age calculation and textual interpretation of the results feature
- Configurable printout format
- Powered directly by the PC / Laptop and automatically launches when connected
- Supplied complete with all accessories in a sturdy carrying case

7154 SpiroUSB with Spirometry PC Software £1024.50

Mouthpieces

7165 Adult Disposable, pk/500 £25.00

7170 Adult Disposable One Way Valve, pk/250 £35.00


Spirometry PC Software (SPCS)



CareFusion



systemone
ONEPATIENTONERECORD



- Compatible with EMIS, SystemOne
- Clear graphic user interface
- Clearer screen layout of results and graphs
- Works on all individual patient basis – so the screens you see are specific to the patient being tested
- Patient specific choice of predicted values
- Results tables can be formatted to appear in the choice of the user
- Predicted values, % predicted and z scores displayed
- Powerful search facilities which allows filtering on patient or groups
- Batch printing options between specific dates
- Increased export functions – to Pdf, MS Word, Excel, XML
- Minimum PC requirement is Windows 2000, Windows XP, Windows Vista
- There is no annual license fee for linking to SystemOne but EMIS users will need to apply for the link with the SPCS software for an annual fee of £150.00

7012 Spirometry PC Software £406.50

One-way Safety Mouthpieces



- Stops patient inhaling through mouthpiece during spirometry test
- Recommend in place of standard mouthpiece as reduces risk of cross infection between patients
- Can be used with most makes of spirometer using 30mm dia mouthpieces

7170 Adult, pk/250 £35.00

MicroGard Pulmonary Filter



- 99.9% bacterial/viral filter
- Designed to stop moisture droplets passing to and from the spirometer
- Used for an inspiratory manoeuvre is performed or when testing patients who may pose a health risk to others
- Can be used with most makes of spirometer (standard 30mm diameter)

7171 Pack/25 £22.95

7172 Pack/100 £79.95

3 Litre Calibration Syringe



Micro Medical spirometers are unaffected by changes in environmental conditions such as temperature, air pressure and humidity, or by moisture build up from the patient's breath during repeated blows.

They can, however, be quickly verified when required by using the 3 litre accuracy check syringe.

7164 Syringe £303.50

Pera® Safe Sterilising Agent



- 10 minute medical device disinfectant/sterilant
- Ideal for effective sterilisation of Micro Medical Gold Standard turbines
- Use 1L of water with 16.2g sachet
- Use 5 litres of water with 81g pot

2386 16.2g Sachet £1.95

2387 81g Pot £7.95

Micro CO



- The gold standard in CO breath testing, and is the most advanced system currently available today
- Features include, fast response time, 1ppm resolution, immediate display of CO levels in PPM and %COHB, colour light indicators
- Free download software
- Supplied complete with a hard shell carry case, sample cardboard mouthpieces, one way plastic mouthpiece adaptor, calibration adaptor, calibration tool, 9V PP3 battery and operating manual
- Dimensions: 17H x 6W x 2.6D cm
- Weight: 160gm
- **Calibration:** In order to keep your Micro CO in pristine working condition, it will require calibration at six monthly intervals. Calibration can easily and safely be carried out by the user with our Medican, 1 litre size calibration cylinder

7146	Micro CO	£245.00
7167	Micro CO & Smoke Check Mouthpieces, pk/250	£19.00
7499	Medican Calibration Kit	£69.95

Baby CO



- Measuring alveolar Carbon Monoxide in PPM concentrations. This is measured simply from the mother's expired breath. From this measurement the amount of Carbon Monoxide in the foetus is derived and is represented as foetal Carboxyhaemoglobin (%FCOHB)
- The Baby CO is the gold standard in CO breath testing and is the most advanced system currently available today. A single breath from the mother will immediately show the percentage of Carbon Monoxide attached to the baby's red blood cells (%FCOHB)
- The mother's CO concentration in PPM will then be displayed, along with the foetal %COHB
- Supplied complete with a hard shell carry case, sample cardboard mouthpieces, one-way plastic mouthpiece adaptor, calibration adaptor, calibration tool, 9V PP3 battery and operating manual
- Free download software
- Dimensions: 17H x 6W x 2.6D cm
- Weight: 160g
- **Calibration:** In order to keep your Baby CO in pristine working condition, it will require calibration at six monthly intervals. Calibration can easily and safely be carried out by the user with our Medican, 1 litre size calibration cylinder

7148	Baby CO	£295.00
7167	MicroCO & Smoke Check Mouthpieces, pk/250	£19.00
7499	Medican Calibration Kit	£69.95

PulmoLife COPD Screening Device



Designed for quick and easy screening of adult smokers for early detection of COPD

- Measures FEV1 — the measurement of choice in COPD screening
- Gives FEV1 result as a percentage of predicted values
- 'Lung Age' Interpretation (Morris and Temple)
- Used as a motivational tool for smoking cessation
- Dimensions: 13.1H x 5.9W x D3.8 cm
- Weight: 96g

7150	PulmoLife	£97.85
------	-----------	--------

Mouthpieces

7165	Adult Disposable, pk/500	£25.00
7170	Adult Disposable One Way Valve, pk/250	£35.00

Smoke Check



- The Smoke Check is the most cost effective CO monitor available today
- Designed as a simple screening test for cigarette consumption
- Provides an instant indication of CO levels in PPM, and backed up with colour light indicators
- Conversion of PPM results to %COHB is easily done using the Smoke Check's smoking cessation guide chart
- Supplied complete with a hard shell carry case, sample cardboard mouthpieces, one way plastic mouthpiece adaptor, calibration adaptor, calibration tool, 9V PP3 battery and operating manual
- Dimensions: 17H x 6W x 2.6D cm
- Weight: 130g
- **Calibration:** In order to keep your Smoke Check in pristine working condition, it will require calibration at six monthly intervals. Calibration can easily and safely be carried out by the user with our Medican, 1 litre size calibration cylinder

7147	Smoke Check	£145.00
7167	MicroCO & Smoke Check Mouthpieces, pk/250	£19.00
7499	Medican Calibration Kit	£69.95



Vitalograph Alpha Touch

NEW

Demo Available



Vitalograph®




- An easy to use, lightweight desktop spirometer that offers the benefits of portability, data integration and built in printer
- Features a Fleisch Pneumotachograph
- Simple icon driven menu with large, high resolution colour touch-screen
- Storage capacity for 10,000 subjects
- Automatic two-way synchronisation of patient data and test sessions between Alpha Touch and Spirotrac
- Over 50 test parameters to choose from
- F/V & V/T curves viewable in real time
- Pre-post bronchodilator comparison
- Single or multi breath
- % predicted or SDS scores comparison
- Hours of battery life for a full day of off-site testing
- Choice of child incentives with sound effects
- Customisable reports including interpretation option, BMI and lung age information
- Supplied with robust carrying case
- 5 Year Warranty subject to annual service contract with Vitalograph



systemone
ONE PATIENT ONE RECORD

EMIS and SystemOne connectivity via Vitalink



emis
accredited partner

9845	Vitalograph Alpha Touch	£1195.00
9846	Vitalograph Alpha Touch with Spirotrac Software	£1495.00
4540	Vitalink® EMIS LV and System One Connectivity Software – annual license fee	£150.00
4541	SafeTway® Mouthpieces, pk/200	£27.95
4543	Bacterial Viral Filters, pk/50	£47.00
4542	Disposable Noseclips, pk/200	£57.00
4545	3 litre Precision Syringe	£325.00
4544	Thermal Paper, pk/50	£72.00

Vitalograph Pneumotrac



PC Not Included



The Vitalograph Pneumotrac is an excellent spirometry choice when long-term storage and analysis of patient results are essential, offering over 30 configurable parameters, storing and trending of test results and export of data into any spreadsheet, database or word processor application.

The Vitalograph Pneumotrac features sophisticated Spirotrac® software and the unparalleled accuracy of Fleisch pneumotachograph flow measuring technology.

- Simply plugs into a PC's USB port for a low-cost spirometry solution
- Adaptable, interactive, user friendly software for full automation
- Network ready, standardised communications with other databases
- Proven spirometry software with automatic test quality control and trending
- Accurate, robust and linear Fleisch pneumotachograph
- Optional Vitalink® Practice Management Connectivity System provides the busy general practice with an automatic two-way interface between the Pneumotrac Spirometry System and EMIS or SystemOne
- 5 Year Warranty subject to annual service contract with Vitalograph

systemone
ONE PATIENT ONE RECORD

EMIS and SystemOne connectivity via Vitalink

Demo Available

emis
accredited partner

4539	Pneumotrac Spirometer	£1095.00
4540	Vitalink® EMIS LV and SystemOne Connectivity Software - annual license fee	£150.00
4541	SafeTWay® Mouthpieces, pk/200	£27.95
4543	Bacterial Viral Filters, pk/50	£47.00
4545	3 litre Precision Syringe	£325.00

Vitalograph In2itive Spirometer



Demo Available

- Combines the benefits of full spirometry, data integration and mobility in a small lightweight handheld spirometer
- Cradle for battery charging, connection to a PC and automatic data exchange saving time and eliminating errors
- Assured hygiene with removable flowhead, SafeTway mouthpieces or BVF Bacterial/Viral Filters
- Hours of battery life for a full day of off-site testing
- Storage capacity for 10,000 subjects
- Over 50 test parameters
- Pre-post bronchodilator comparison
- Single and multi-breath
- F/V or V/T curves viewable in real time
- Choice of child incentives with sound effects
- Icon driven menus
- Print via your PC or via Spirotrac or print direct to a USB A4 letter size printer using optional printer cradle
- Supplied with range of accessories and robust carrying case
- 5 Year Warranty subject to annual service contract with Vitalograph



8938	in2itive Lite with Vitalograph Reports Software	£1095.00
8937	in2itive with Printer Cradle and A4 Printer	£1295.00
8936	in2itive with PC Cradle and Spirotrac V Software	£1295.00
4540	Vitalink® EMIS LV and Systeem One Connectivity Software – annual license fee	£150.00
4541	SafeTway® Mouthpieces, pk/200	£27.95
4543	Bacterial Viral Filters, pk/50	£47.00
4542	Disposable Noseclips, pk/200	£57.00
4545	3 litre Precision Syringe	£325.00



Vitalograph Alpha Spirometer

The Vitalograph ALPHA provides a wide range of measurements, including VC, FVC, FEV1, FEV1%, PEF and FEF 25-75. Information is presented in an easy-to-read format on a clear, bright colour screen and multiple pre-tests can be stored for subsequent post-testing. Featuring the unparalleled accuracy of a Fleisch pneumotachograph and linked to your PC via USB connection, the ALPHA offers a range of additional features including integral battery and mains power options, report generation software, real time clock, easy to clean flowhead, smart carry case and a fast, quiet printer that prints 140 tests per roll.

Demo Available

- Dedicated all-in-one desktop unit
- Simple to use
- Rapid and low cost testing
- Clear and bright colour screen
- Accurate Fleisch pneumotachograph
- Measures VC, FVC, FEV1, FEV1%, PEF, FEF 25-75 and more
- 5 Year Warranty subject to annual service contract with Vitalograph



1

4538	Alpha Spirometer	1	£1095.00
4541	SafeTWay® Mouthpieces, pk/200		£27.95
4543	Bacterial Viral Filters, pk/50		£47.00
4544	Thermal Paper, pk/5		£72.00
4545	3 litre Precision Syringe	2	£325.00



2



Vitalograph asma-1 & asma-1 usb Electronic Asthma Monitors



The Vitalograph asma-1 is an easy to use, yet sophisticated, electronic asthma monitor that makes the monitoring of asthma simpler and more accurate than ever before. By measuring FEV1, as well as peak flow, the asma-1 gives a clearer picture of lung health. Storage capacity of up to 600 test sessions and automatically assesses test quality.

- Choice of two monitors, both with identical specification, the model asma-1usb links to your PC to produce pdf and hard copy colour reports
- Measures FEV1 and peak flow
- Simple to use
- Large, easy to read display
- No need for record cards
- Quality of blow indicator
- Suitable for multi-subject use with standard size SafeTWay® mouthpieces
- Uses long life AAA batteries

**AWARD
WINNING
DESIGN
PRODUCT**

4546	asma-1	1	£35.00
4547	asma-1usb	2	£135.00
4541	SafeTWay® Mouthpieces, pk/200		£27.95
4542	Disposable Noseclips, pk/200		£57.00

Vitalograph copd-6 & copd-6 usb COPD Screening Devices



The Vitalograph copd-6 is an accurate and cost effective 'case selection' or screening device that identifies potential cases of COPD

- Choice of two monitors, both with identical specification, the model copd-6usb links to your PC to produce pdf and hard copy colour reports
- Identifies those at risk of COPD at the pre-symptomatic stage
- Screens out normal FEV1, without the risk of false negatives
- Displays the key parameters for clinical interpretation including obstructive index and COPD classification (stage I - IV)
- Displays FEV1, FEV6, ratio and % predicted, obstructive index, COPD classification and lung age
- Detachable flowhead makes cleaning simple
- Built-in quality of blow indicator
- Designed for multi-subject use

**AWARD
WINNING
DESIGN
PRODUCT**

4548	copd-6	3	£69.50
4549	copd-6usb	4	£135.00
4541	SafeTWay® Mouthpieces, pk/200		£27.95
4542	Disposable Noseclips, pk/200		£57.00

Vitalograph Peak Flow Meter

**Unbeatable
Value**



- Accurate
- Compact, lightweight and portable
- Robust and hard wearing
- Features comfortable integral mouthpiece
- Easily read scale markers
- Suitable for adults and children
- Suitable for multi-subject use with SafeTWay® mouthpiece

4550	Peak Flow Meter	£4.95
4541	SafeTWay® Mouthpieces, pk/200	£27.95

Vitalograph Consumables

Vitalograph SafeTWay® Mouthpieces

- For use with peak flow meters and spirometers for the effective control of cross infection risk
- One-way valve
- Disposable
- A low cost alternative to BVF for expiratory testing



4541	SafeTWay® Mouthpieces, pk/200	£27.95
------	-------------------------------	--------

Vitalograph Disposable Noseclips

- Disposable nose clips for spirometry testing
- Designed for optimum comfort



4542	Disposable Noseclips, pk/200	£57.00
------	------------------------------	--------

Clement Clarke Mini-Wright Peak Flow



- Choice of two models; Standard or Low Range
- Standard: Range 60-800 l/min
- Low Range: 30-400 l/min, suitable for children and patients with severe respiratory conditions
- Both devices conform to European Standards
- Supplied with instructions which include patient recording chart

4950	Mini-Wright Standard EU	£8.25
4951	Mini-Wright EU Low Range	£8.25
4969	Adult Disposable One Way Mouthpieces Expiratory, pk/200	£29.75
4971	Paediatric Mouthpiece Adaptor for One Way Mouthpieces, pk/20	£2.00
4967	Adult Disposable Two Way Mouthpieces, pk/500	£12.50
4968	Paediatric Disposable Two Way Mouthpieces, pk/500	£19.95

Clement Clarke One Flow FVC Memo Spirometer



Measures PEF, FEV1, FVC and FEV1/FVC Ratios

- Ideal for patient screening in hospitals, GP surgeries, at home or for occupational health
- One Flow FVC Kit includes PC software, USB cable and a spare mouthpiece
- Software enables flow/volume and volume/time curves to be produced with patient data in tabular format
- Reversibility testing can be performed with 3 pre and 3 post measurements and superimposed curves
- Memory Storage: 100 PEF and 100 FEV1
- Easy to use, portable, reliable and accurate
- Supplied with 4 x AAA batteries

4965	One Flow FVC Memo Spirometer	1	£250.00
4966	One Flow FVC Memo Kit Spirometer	2	£300.00
4969	Adult Disposable One Way Mouthpieces Expiratory, pk/200		£29.75
4967	Adult Disposable Two Way Mouthpieces, pk/500		£12.50



Clement Clarke Mini-Wright Digital Peak Flow

- Mini-Wright Digital not only measures Peak Expiratory Flow (PEF) but also measures FEV₁ for both Asthma and COPD patients
- Designed to monitor asthma
- Screening device for COPD patients
- Measures both Peak Expiratory Flow (PEF) and Forced Expiratory Volume (FEV₁)
- Multi patient use with One Way Valve Mouthpieces
- FEV₁ and PEF zoning can be set
- Capacity to store 240 PEF and 240 FEV₁ results
- Water tight seal making MDD waterproof for simple sterilisation
- Just one button is needed to use device, making it simple to use
- Ideal device to use as electronic Peak Flow diary
- Compact, lightweight

Mini Wright Digital Software

- Allows viewing and trending results for effective monitoring of Asthma and COPD
- Data downloaded via USB port and Mini-Wright Digital Download cradle
- Date and Time stamping eliminates the need for record cards
- Viewed parameters are available for tables and graphs
- Graphical data selectable morning, afternoon or day
- MDD software configure Mini-Wright digital with zones and reference values
- Allows a build up of patient database



4952	Mini-Wright Digital Download with Cradle & Software	£75.00
4969	Adult Disposable One Way Mouthpieces Expiratory, pk/200	£29.75
4967	Adult Disposable Two Way Mouthpieces, pk/500	£12.50

Clement Clarke Mouthpieces

- Mouthpieces are colour-coded
- White with Blue lettering for inhalation only
- White with Red lettering for exhalation only
- Large arrows on the outside of the tube identifying which way air can flow

4969	Adult Disposable One Way Expiratory Red, pk/200	£29.75
4970	Adult Disposable One Way Inspiratory Blue, pk/200	£35.00
4971	Paediatric Mouthpiece Adaptor for One Way, pk/20	£2.00
4967	Adult Disposable Two Way, pk/500	£12.50
4968	Paediatric Disposable Two Way, pk/500	£19.95



Nebulair Hi-Flow Nebuliser



- High output mains nebuliser
- Maintenance free piston compressor
- Compact unit with built in carry handle
- Nebulises all modern medications efficiently
- Includes 2 sets of accessories, with child and adult masks, chambers, filters, tubing and mouthpieces
- Can be used lying down
- Output 3.5 bar
- Weight: 2.4kg
- Very quiet operation



6936	Nebulair Hi-Flow Nebuliser	£85.50
6932	Replacement Adult Face Mask, Chamber, Tubing Kit	£2.75
6933	Replacement Child Face Mask, Chamber, Tubing Kit	£2.75
6934	Replacement Universal Mouthpiece Nebuliser Kit	£2.75

Guardian F4000 Ned-Aid



- Modern high efficiency compact nebuliser
- Unobtrusive and quiet piston compressor
- Black carry and storage case with handle and strap
- Max pressure: 1.5 bar approx, Compressor output: 9 litres per minute
- Dimensions: 18H x 14W x 13D cm
- Weight: 1.3kg

Supplied with: Adult soft face mask, Child soft face mask, Connecting tubing, Nebulising chamber, Mouthpiece and nasal adaptor

6935	F4000 Ned-Aid	£57.50
6932	Replacement Adult Face Mask, Chamber, Tubing Kit	£2.75
6933	Replacement Child Face Mask, Chamber, Tubing Kit	£2.75
6934	Replacement Universal Mouthpiece Nebuliser Kit	£2.75

Lella-the-Ladybug Childs Nebuliser



- Perfect For Small Children
- Quiet Operation
- Zip carry case is transparent and shaped to fit Lella and all her accessories
- 'Lella the Ladybug' Story Book
- Suitable for Bronchodilators - And other medications held in suspension
- Medication capacity 2ml (min) - 8ml (max)
- Operating measure: 0.9bar
- Dimensions 21H x 16W x 13D cm
- Weight: 1.6kg

Supplied with: Story Book, Rhino Clear Nasal Washer, Chamber, Large and small face masks, Mouthpiece, Tubing, Carry case

6937	Lella-the-Ladybug Childs Nebuliser	£89.00
------	------------------------------------	--------

Omron Comp Air C28 Nebuliser



- Compact design
- Wide medication compliance
- Simple one-button operation
- Quick nebulisation times
- Handy carrying case
- Nebuliser rate average 0.4ml/minute
- Dimensions: 10.3H x 17WD x 18.2D cm
- Weight: 1.9kg

4857	Omron Comp Air C28 Nebuliser	£57.50
------	------------------------------	--------

Nebuliser Set



- Patients mask with nebuliser chamber and tubing
- Nebuliser mouthpiece set for those patients who use mouthpiece instead of face mask

6932	Nebuliser Set Adult	1	£2.75
6933	Nebuliser Set Child		£2.75
6934	Nebuliser Set Mouthpiece, Adult	2	£2.75

Medix Actineb



- Suitable for patients of all ages
- Fast and effective therapy for bronchodilators (relievers), corticosteroids (preventers) and antibiotics
- Nebulisation Time: 4.5 to 7.5 minutes for 2.5ml
- Discreet storage compartment allows for all accessories and cable to be neatly stored
- Supplied with Nebuliser chamber, driveline, inlet filter, adult face mask
- Dimensions: 28H x 20W x 12.5D cm
- Weight: 2.3kg

4941	Actineb	£64.00
4946	Actineb Adult Year Pack	£20.00
4947	Actineb Child Year Pack	£20.00

Medix AC2000



- Fast and effective therapy for bronchodilators (relievers), corticosteroids (preventers) and antibiotics
- Easy to use & suitable for patients of all ages
- Nebulisation Time: 2.5 - 5.5 minutes for 2.5ml
- Continuous self cooling motor
- Discreet storage compartment allows for all accessories and cable to be neatly stored
- Supplied with Nebuliser chamber, driveline, inlet filter, adult face mask
- Dimensions: 37.5H x 23.5W x 11.5D cm
- Weight: 2.45kg

7582	AC2000	£114.00
4944	AC2000/World Traveller Adult Year Pack	£20.00
4945	AC2000/World Traveller Child Year Pack	£20.00

Medix Microneb



Handheld, portable, silent nebuliser - ideal product for those who need to nebulise on the move

- Ideal for children and adults
- Can be used by Asthma and COPD sufferers
- Ensures low particle size which improves medication delivery
- Simple one button operation
- Battery powered, 2 AA batteries can last for up to 8 days - can be used with rechargeable batteries
- Optional mains adaptor for use in home
- Lightweight - Weighs just 300g
- Year Pack contains: 4 x adult or 4 x child face masks and 1 x mouthpiece

4937	Microneb	£110.00
4938	Microneb Mains Adaptor	£39.00
4939	Microneb Adult Year Pack	£10.00
4940	Microneb Child Year Pack	£10.00

Medix Econoneb



- Designed for continuous use making it a reliable nebuliser for Clinic, Hospital and Home use
- Suitable for patients of all ages
- Nebulises all commonly used respiratory drugs
- Easy to use controls and carrying handle
- High flow rate which is suitable for nebulising bronchodilators (relievers), corticosteroids (preventers) and antibiotics
- Supplied with Nebuliser chamber, driveline, inlet filter, adult face mask
- Dimensions: 21H x 18.5W x 18.5D cm
- Weight: 2.35kg

9436	Econoneb	£99.00
9437	Econoneb Adult Year Pack	£20.00
9438	Econoneb Child Year Pack	£20.00

Medix Nebuliser Year Packs



- Year pack contains: 4 x nebuliser chambers, 4 x drivelines, 4 x adult or 4 x child face masks, 1 x mouthpiece and 4 x inlet discs

4944	AC2000/World Traveller Adult Year Pack	£20.00
4945	AC2000/World Traveller Child Year Pack	£20.00
4946	Actineb Adult Year Pack	£20.00
4947	Actineb Child Year Pack	£20.00
4939	Microneb Adult Year Pack	£10.00
4940	Microneb Child Year Pack	£10.00
9437	Econoneb Adult Year Pack	£20.00
9438	Econoneb Child Year Pack	£20.00

Medix World Traveller

The World Traveller operates from any mains supply worldwide. It can also be used in the car with a 12V DC lead



- Fast and effective therapy every time for bronchodilators (relievers), corticosteroids (preventers) and antibiotics
- Nebulisation Time: 2.5 - 5.5 minutes for 2.5ml
- Suitable for patients of all ages
- Optional battery pack available
- Easy carrying handle
- Supplied with mains & 12v leads, Nebuliser chamber, driveline, inlet filter, adult face mask
- Dimensions: 37.5H x 23.5W x 11.5D cm
- Weight: 2.5kg (without battery) 2.8kg (with battery)

4943	World Traveller	£240.00
5568	World Traveller Rechargeable Battery Pack	£60.00
4944	AC2000/World Traveller Adult Year Pack	£20.00
4945	AC2000/World Traveller Child Year Pack	£20.00

